


LC|GC

north america

solutions for separation scientists

Volume 32 Number 3 March 2014



SUPELCO[®]
Solutions within.™

NEW! Titan™ UHPLC Columns Starting under \$300

Titan UHPLC columns provide the performance scientists expect (typically >250,000 N/m), but at a greatly reduced cost. Get yours today!

TITAN™

Totally Porous, Monodisperse Silica

Explore the solutions within



For new product information,
ordering and real time availability, visit
sigma-aldrich.com/titan

More Details >

SIGMA-ALDRICH[®]

Outperforms other UHPLC columns with the narrowest particle size distribution of any sub-2 µm totally porous particles

The Titan™ C18 UHPLC columns are based on 1.9 µm totally porous, monodisperse silica particles. These particles are the result of a newly developed, patent-pending, Ecoporous™ silica manufacturing process.

Titan delivers:

- Efficiency > 250,000 N/m, comparable to leading UHPLC columns
- Narrow particle size distribution translates to rugged, reliable performance
- A great value, available at half the cost of other UHPLC columns



For new product information, ordering and real time availability, visit sigma-aldrich.com/titan

Titan C18 Performance Comparison in Acetonitrile

Comparison of the Titan performance to commercially available sub-2 µm columns


Dionex® 3000 (Low D tubing)

column: 5 cm x 3.0 mm I.D.
mobile phase: 60% acetonitrile
temp: 35 °C
flow rate: 0.9 mL/min (4 mm/s)
detection: 254 nm

1. Uracil
2. Diazepam
3. Toluene
4. Naphthalene
5. Biphenyl


Titan C18, 1.9 µm

N= 14,710
N/m= 294,200
h= 1.75
σ= 1.05
Pressure*= 4,100psi




Competitor A C18, 1.8 µm

N= 9,260
N/m= 185,200
h= 3.18
σ= 1.02
Pressure*= 4,650psi



Competitor B C18, 1.7 µm

N= 9,783
N/m= 195,700
h= 2.84
σ= 1.16
Pressure*= 4,900psi



* Total System Pressure

LCGC

north america

solutions for separation scientists

Volume 32 Number 3 March 2014
www.chromatographyonline.com



Analyzing Amphetamines and Bath Salts with SPE and LC-MS-MS

In-Capillary Sample Concentration in CE
LC Method Scaling with Gradient Separations
Safe Storage and Use of GC Gases

THOUSANDS
OF PARTS & SUPPLIES

HUNDREDS
OF INSTRUMENTS SERVICED

FIRST
NAME TO CALL

Agilent CrossLab delivers quality parts and supplies, tested to ensure seamless performance with major brands of GC and HPLC instruments. Our expert service engineers are certified to repair, maintain and ensure regulatory compliance of instrument platforms from all major manufacturers. And with the CrossLab Stand Behind Warranty, you get confidence, not compromise.

The Measure of Confidence

© Agilent Technologies, Inc. 2014



**STAND
BEHIND
WARRANTY**

Learn why we should be
your trusted source at
agilent.com/chem/CrossLab

Agilent
CrossLab
SUPPLIES | SERVICES



Agilent Technologies



HPLC Polymers 101

The power of polymeric columns

Polymer HPLC columns have a lot of benefits. They don't require any functionalization for reversed-phase separations, and rigid polymeric supports intrinsically resist chemical and pH degradation, a fundamental problem with silica columns. Plus, polymer's inertness to most chemical environments makes it a robust and economical solution.

Hamilton offers a line of pH stable polymer HPLC columns for reversed phase, anion exchange, cation exchange and ion exclusion separations perfect for pharmaceuticals, small molecules, proteins, peptides, DNA, organic and inorganic ions and more.

Advantages of polymer over silica

- ✓ pH range of 1–13
- ✓ Widest chemical compatibility
- ✓ Temperatures higher than 60 °C
- ✓ Maximum sample recovery
- ✓ Longest average life span

To learn more about how polymer columns can perform for you, visit www.ham-info.com/0805-1 or call toll free **1-888-525-2123**.

© 2014 Hamilton Company. All rights reserved.
Images Copyright Rangizzz and Carolina K. Smith, M.D., 2014
Used under license from Shutterstock.com



Chem Service
inc.

Environmental & Analytical Standards



ISO Guide 34:2009 Accredited
ISO 17025:2005 Accredited
ISO 9001:2008 Certified

- Over 13,000 Reference Standards
- Custom Synthesis
- Quick Delivery

STANDARDS YOU CAN RELY ON!

- **High Quality:** High purity chemicals for use as certified reference materials.
- **Cost Effective:** Products packaged in small quantities to minimize expenses.
- **Wide Selection:** 1,200+ pesticide and metabolite reference standards (including banned and discontinued items).
- **Certification backed with experience:** Over 50 years of experience. ISO Guide 34:2009, ISO 9001:2008, and ISO/IEC 17025:2005 certification.

Toll Free: 1-800-452-9994
Direct: 610-692-3026 • Fax: 610-692-8729



ORDER ONLINE
www.chemservice.com

New CATALOG!
request a FREE copy
(Print, CD-ROM or electronic version)



LCGC[®]

north america

MANUSCRIPTS: For manuscript preparation guidelines, see chromatographyonline.com/AuthorInfo, or call The Editor, (732) 596-0276. LCGC welcomes unsolicited articles, manuscripts, photographs, illustrations, and other materials but cannot be held responsible for their safekeeping or return. Every precaution is taken to ensure accuracy, but LCGC cannot accept responsibility for the accuracy of information supplied herein or for any opinion expressed.

SUBSCRIPTIONS: For subscription and circulation information: LCGC, P.O. Box 6168, Duluth, MN 55806-6168, or call (888) 527-7008 (7:00 a.m.–6:00 p.m. central time). International customers should call +1-218-740-6477. Delivery of LCGC outside the United States is 14 days after printing. For single and back issues, call (800) 598-6008 or (218) 740-6480. (LCGC Europe and LCGC Asia Pacific are available free of charge to users and specifiers of chromatographic equipment in Western Europe and Asia and Australia, respectively.)

CHANGE OF ADDRESS: Send change of address to LCGC, P.O. Box 6168, Duluth, MN 55806-6168; alternately, send change via e-mail to magazines@superfill.com or go to the following URLs:

• Print: <https://advanstar.replycentral.com/Default.aspx?PID=469>

• Digital: <https://advanstar.replycentral.com/?PID=469&V=DIGI>

Allow four to six weeks for change. PUBLICATIONS MAIL AGREEMENT No. 40612608. Return all undeliverable Canadian addresses to: IMEX Global Solutions, P.O. Box 25542, London, ON, N6C 6B2, CANADA. Canadian GST number: R-124213133RT001.

DIRECT MAIL LIST RENTAL: Contact Tamara Phillips, tel. (440) 891-2773, e-mail TPhillips@advanstar.com.

REPRINTS: Reprints of all articles in this issue and past issues of this publication are available (500 minimum). Call 877-652-5295 ext. 121 or e-mail bkolb@wrights-media.com. Outside US, UK, direct dial: 281-419-5725. Ext. 121

MARKETING DEVELOPMENT/CLASSIFIED: Contact Tod McCloskey, tel. (440) 891-2739, fax (440) 826-2865.

RECRUITMENT: Contact Tod McCloskey, tel. (440) 891-2739, fax (440) 826-2865.

INTERNATIONAL LICENSING: Contact Maureen Cannon, tel. (440) 891-2742, fax (440) 891-2650, or e-mail mcannon@advanstar.com.

selective GC detection for complex matrices



CATALYTIC COMBUSTION IONIZATION (TID-10 Ceramic Ion Source & Oxidizing Detector Gas)

CCID - transient flame ignition as compounds with chains of CH₂ groups impact a hot TID-10 ceramic surface - selective for Alkanes, Alkenes, FAMES, Triglycerides, etc., but no Aromatic or Cyclo-HC response

selectivity for Saturate vs. Unsaturate Carbon bonds is achieved by adjusting Oxygen level in the detector gas

surface ionization on TID-10 ceramic also adds responses for Oxygenates & Pyrrole functional groups.

GASOLINE ANALYSIS



**DETECTOR
Engineering &
Technology, inc.**

486 N. Wiget Lane, Walnut Creek, CA 94598
tel: 925-937-4203, email: detplp@aol.com

www.det-gc.com

BPA
WORLDWIDE[™]
BUSINESS

**american
business
media**
The Association of Business Information Companies

©2014 Advanstar Communications Inc. All rights reserved. No part of this publication may be reproduced or transmitted in any form or by any means, electronic or mechanical including by photocopy, recording, or information storage and retrieval without permission in writing from the publisher. Authorization to photocopy items for internal/educational or personal use, or the internal/educational or personal use of specific clients is granted by Advanstar Communications Inc. for libraries and other users registered with the Copyright Clearance Center, 222 Rosewood Dr. Danvers, MA 01923, 978-750-8400 fax 978-646-8700 or visit <http://www.copyright.com> online. For uses beyond those listed above, please direct your written request to Permission Dept. fax 440-756-5255 or email: mcannon@advanstar.com.

Advanstar Communications Inc. provides certain customer contact data (such as customer's name, addresses, phone numbers, and e-mail addresses) to third parties who wish to promote relevant products, services, and other opportunities that may be of interest to you. If you do not want Advanstar Communications Inc. to make your contact information available to third parties for marketing purposes, simply call toll-free 866-529-2922 between the hours of 7:30 a.m. and 5 p.m. CST and a customer service representative will assist you in removing your name from Advanstar's lists. Outside the U.S., please phone 218-740-6477.

LCGC North America does not verify any claims or other information appearing in any of the advertisements contained in the publication, and cannot take responsibility for any losses or other damages incurred by readers in reliance of such content.

To subscribe, call toll-free 888-527-7008. Outside the U.S. call 218-740-6477.

Advanstar Communications Inc. (www.advanstar.com) is a leading worldwide media company providing integrated marketing solutions for the Fashion, Life Sciences and Powersports industries. Advanstar serves business professionals and consumers in these industries with its portfolio of 91 events, 67 publications and directories, 150 electronic publications and Web sites, as well as educational and direct marketing products and services. Market leading brands and a commitment to delivering innovative, quality products and services enables Advanstar to "Connect Our Customers With Theirs." Advanstar has approximately 1000 employees and currently operates from multiple offices in North America and Europe.



Bruker Chemical Analysis

● GC, GC-MS, LC-MS & ICP-MS



SCION GC GC Systems

SCION SQ & TQ GC-MS Systems

aurora Elite ICP-MS Systems

EVOQ LC-MS LC-MS Systems

Innovative Chromatography and Mass Spectrometry Systems for Chemical Analyses

- Robust, powerful and easy-to-implement analytical system solutions for any challenge
- Integrated instrument and software systems optimized to enhance the productivity and capabilities of any analytical laboratory
- Market leading portfolio of technologies, tools and solutions including GC, GC-Single Quads, GC-Triple Quads, LC-Triple Quads, ICP-MS, ESI-Q-TOF, LC-IonTrap, FTMS and MALDI-TOF/TOF
- Consistently and expertly supported to ensure you meet all your goals



Contact us for more details at cam-sales@bruker.com


Innovation with Integrity

GC, GC-MS, LC-MS, ICP-MS

LCGC®

north america

solutions for separation scientists



Cover photography by Joe Zucgic,
Joe Zucgic Photography

Cover materials courtesy of Agilent
Technologies, Imtakt USA, and
UCT, LLC

WEB SEMINARS

Increasing Laboratory Efficiency through a Harmonized Compliance Solution

Gary Powers and Shae Taylor,
Agilent Technologies

chromatographyonline.com/WebSeminar

NEW ON LCGCtv

Koen Sandra
on analyzing biopharmaceuticals
with LC-MS

Dwight Stoll
on 2D LC

Mary Ellen McNally
on the proper use of an
internal standard

Michael Holcapek
on trends in LC-MS

Frederic Lynen
on analyzing drugs and metabolites
in biological matrices

www.learnpharmascience.com/lcgc

Like LCGC on Facebook:
www.facebook.com/lcgc magazine

Follow LCGC on Twitter:
https://twitter.com/LC_GC

Join the LCGC Group on LinkedIn
<http://linkd.in/LCGCgroup>



CONTENTS

COLUMNS

174 COLUMN WATCH

In-Capillary Sample Concentration in CE

Michael C. Breadmore and Cari E. Sanger-van de Griend
Stacking and sweeping injection methods are compared
and practical suggestions are provided about how to
enhance sensitivity for charged and neutral compounds.

188 LC TROUBLESHOOTING

LC Method Scaling, Part II: Gradient Separations

John W. Dolan

If scaling isocratic separations is so simple, why is gradient
scaling so confusing?

194 GC CONNECTIONS

The Storage and Use of Gases for Gas Chromatography

John V. Hinshaw

A discussion on the safe storage and use of calibration,
carrier, and detector gas cylinders for small or large organizations

226 THE ESSENTIALS

Mixed-Mode HPLC Separations: What, Why, and How

The mechanisms of mixed-mode chromatography are
explained, with reference to the parameters that are
used to optimize selectivity and retention.

PEER-REVIEWED ARTICLES

200 Analysis of Amphetamines and Synthetic Cathinones in Hair Samples Using LC-Tandem Mass Spectrometry and Solid-Phase Extraction

Albert A. Elian, Jeffery Hackett, and Kerrie T. Donovan

An efficient method for forensic analysis of amphetamines and
synthetic cathinones — the illicit drugs often called “bathsalts” —
in hair samples

210 Application of Pyrolysis-Gas Chromatography-Mass Spectrometry for the Identification of Polymeric Materials

*Peter Kusch, Gerd Knupp, Wolfgang Fink, Dorothee Schroeder-Obst,
Volker Obst, and Johannes Steinhaus*

This pyrolysis-GC-MS method enables direct analysis of solid or
liquid polymers without sample pretreatment, as illustrated here for
various materials, including a dental filling material and a car
wrapping foil.

DEPARTMENTS

Peaks of Interest 172 Products & Resources .218 Ad Index 224

LCGC North America (ISSN 1527-5949 print) (ISSN 1939-1889 digital) is published monthly with 1 additional issue in August as Buyers Guide by Advanstar Communications Inc., 131 West First Street, Duluth, MN 55802-2065, and is distributed free of charge to users and specifiers of chromatographic equipment in the United States and Canada. Single copies (prepaid only, including postage and handling): \$15.50 in the United States, \$17.50 in all other countries; back issues: \$23 in the United States, \$27 in all other countries. LCGC is available on a paid subscription basis to nonqualified readers in the United States and its possessions at the rate of: 1 year (13 issues), \$74.95; 2 years (26 issues), \$134.50; in Canada and Mexico: 1 year (13 issues), \$95; 2 years (26 issues), \$150; in all other countries: 1 year (13 issues), \$140; 2 years (26 issues), \$250. Periodicals postage paid at Duluth, MN 55806 and at additional mailing offices. POSTMASTER: Please send address changes to LCGC, P.O. Box 6168, Duluth, MN 55806-6168. PUBLICATIONS MAIL AGREEMENT NO. 40612608, Return Undeliverable. Canadian Addresses to: IMEX Global Solutions, P. O. Box 25542, London, ON N6C 6B2, CANADA Canadian GST number: R-124213133RT001. Printed in the USA.

BRING EVERYTHING IN YOUR LABORATORY TOGETHER.
AND WE MEAN EVERYTHING.



Introducing the world's most comprehensive data management platform.

Your instruments. Your systems. Your samples. Your people. NuGenesis® 8 integrates them all. That's data management, eLab notebook and laboratory execution capabilities on a single platform. Now science-driven, business-minded organizations like yours can capture, catalog, access and deliver information with more insight than ever before. Easy access to this insight means better collaboration. And, in the end, greater results.


To experience the power of NuGenesis 8 for yourself, visit waters.com/nugensis8.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Pharmaceutical & Life Sciences | Food | Environmental | Clinical | Chemical Materials

Leave the others behind, try
something new....
in Purge & Trap



www.cdsanalytical.com

LCGC[®]
north america

Publishing & Sales

485F US Highway One South, Suite 210, Iselin, NJ 08830
tel. (732) 596-0276 fax (732) 647-1235

Science Group Publisher **Michael J. Tessalone**
MTessalone@advanstar.com
Associate Publisher **Edward Fantuzzi**
EFantuzzi@advanstar.com
East Coast Sales Manager **Stephanie Shaffer**
SShaffer@advanstar.com

Editorial

Editorial Director **Laura Bush**
LBush@advanstar.com
Managing Editor **Megan L'Heureux**
MLheureux@advanstar.com
Group Technical Editor **Stephen A. Brown**
SBrown@advanstar.com
Associate Editor **Cindy Delonas**
CDelonas@advanstar.com
Art Director **Dan Ward**
DWard@media.advanstar.com

Vice President Sales **Russell Pratt**
RPratt@advanstar.com

Marketing Manager **Anne Young**
AYoung@advanstar.com

Classified/Recruitment Sales Representative **Tod McCloskey**
TMcCloskey@advanstar.com

Direct List Rental Sales **Tamara Phillips**
TPhillips@advanstar.com

Permissions **Maureen Cannon**
MCannon@advanstar.com

Reprint Services **877-652-5295 ext. 1211** bkolb@wrightsmedia.com
Outside US, UK, direct dial: 281-419-5725. Ext. 121

Production Manager **Jesse Singer**
JSinger@media.advanstar.com

Assistant Audience Development Manager **Gail Mantay**
GMantay@advanstar.com

ADVANSTAR
SCIENCE

Chief Executive Officer **Joe Loggia**

Chief Executive Officer Fashion Group, Executive Vice-President **Tom Florio**

Executive Vice-President, Chief Administrative Officer &
Chief Financial Officer **Tom Ehardt**

Executive Vice-President **Georgiann DeCenzo**

Executive Vice-President **Chris DeMoulin**

Executive Vice-President **Ron Wall**

Executive Vice-President, Business Systems **Rebecca Evangelou**

Executive Vice-President, Human Resources **Julie Molleston**

Sr Vice-President **Tracy Harris**

Vice-President, Information Technology **Joel Horner**

Vice-President, Legal **Michael Bernstein**

Vice-President, Media Operations **Francis Heid**

Vice-President, Treasurer & Controller **Adele Hartwick**

PINNACLE PCX
(Analyze with Power)

Our Pinnacle PCX, with the most trusted
methodology for amino acid analysis, enables
decidedly strong results.

www.pickeringslabs.com

PICKERING
LABORATORIES
CATALYST FOR SUCCESS

Rediscover Your Curiosity



“I wonder what happens if...”

Remember your first discovery? That moment when you first realized you had found something new? Curiosity drives science, and real breakthroughs depend on analytical tools that provide accurate, reliable data. We develop innovative, high-quality chromatography

products and provide expert technical support, so you can trust your results and focus on asking the questions that lead to new advances. With Restek, you can *rediscover your curiosity and push the frontiers of science.*


RESTEK
Rediscover Restek

Want to learn more about how Restek products can help you advance your science?

Sign up to receive free literature in your area at www.restek.com/rediscover



Editorial Advisory Board

Kevin D. Altria GlaxoSmithKline, Ware, United Kingdom
Jared L. Anderson The University of Toledo, Toledo, Ohio
Daniel W. Armstrong University of Texas, Arlington, Texas
Michael P. Balogh Waters Corp., Milford, Massachusetts
Brian A. Bidlingmeyer Agilent Technologies, Wilmington, Delaware
Dennis D. Blevins Agilent Technologies, Wilmington, Delaware
Peter Carr Department of Chemistry, University of Minnesota, Minneapolis, Minnesota
Jean-Pierre Chervet Antec Leyden, Zoeterwoude, The Netherlands
John W. Dolan LC Resources, Walnut Creek, California
Michael W. Dong Genentech, San Francisco, California
Roy Eksteen Sigma-Aldrich/Supelco, Bellefonte, Pennsylvania
Anthony F. Fell School of Pharmacy, University of Bradford, Bradford, United Kingdom
Francesco Gasparrini Dipartimento di Studi di Chimica e Tecnologia delle Sostanze Biologicamente Attive, Università "La Sapienza," Rome, Italy
Joseph L. Glajch Momenta Pharmaceuticals, Cambridge, Massachusetts
Davy Guillaume University of Geneva, University of Lausanne, Geneva, Switzerland
Richard Hartwick PharmAssist Analytical Laboratory, Inc., South New Berlin, New York
Milton T.W. Hearn Center for Bioprocess Technology, Monash University, Clayton, Victoria, Australia
Emily Hilder University of Tasmania, Hobart, Tasmania, Australia
John V. Hinshaw BPL Global, Ltd., Hillsboro, Oregon
Kiyokatsu Jinno School of Materials Science, Toyohashi University of Technology, Toyohashi, Japan

Ira S. Krull Northeastern University, Boston, Massachusetts
Ronald E. Majors LCGC columnist and analytical consultant, West Chester, Pennsylvania
R.D. McDowall McDowall Consulting, Bromley, United Kingdom
Michael D. McGinley Phenomenex, Inc., Torrance, California
Victoria A. McGuffin Department of Chemistry, Michigan State University, East Lansing, Michigan
Mary Ellen McNally E.I. du Pont de Nemours & Co., Wilmington, Delaware
Imre Molnár Molnar Research Institute, Berlin, Germany
Glenn I. Ouchi Brego Research, San Jose, California
Colin Poole Department of Chemistry, Wayne State University, Detroit, Michigan
Fred E. Regnier Department of Chemistry, Purdue University, West Lafayette, Indiana
Pat Sandra Research Institute for Chromatography, Kortrijk, Belgium
Peter Schoenmakers Department of Chemical Engineering, University of Amsterdam, Amsterdam, The Netherlands
Kevin Schug University of Texas, Arlington, Texas
Dwight Stoll Gustavus Adolphus College, St. Peter, Minnesota
Michael E. Swartz Ariad Pharmaceuticals, Cambridge, Massachusetts
Thomas Wheat Waters Corporation, Milford, Massachusetts

CONSULTING EDITORS: Jason Anspach, Phenomenex, Inc.; Stuart Cram, Thermo Fisher Scientific; David Henderson, Trinity College; Tom Jupille, LC Resources; Sam Margolis, The National Institute of Standards and Technology; Joy R. Miskis, Bioanalytical Solutions LLC; Frank Yang, Micro-Tech Scientific.

SUPELCO®
Solutions within.™

Reveal the truth

Improve results, reduce cost and increase throughput

Proprietary HybridSPE®-Phospholipid technology removes proteins and phospholipids from biological samples, reducing the need for costly and time-consuming column replacement while improving LC-MS data accuracy and reproducibility.



For videos and product information, visit
sigma-aldrich.com/hybridspe-pl



©2014 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. HybridSPE is a registered trademark of Sigma-Aldrich Co. LLC. Solutions within is a trademark of Sigma-Aldrich Co. LLC.

12069

SIGMA-ALDRICH®

Experience the quintessence

Milli-Q® Integral system
pure and ultrapure water at your fingertips.

- Dual POD (point-of-delivery) concept saves space and increases convenience.
 - Lower running costs and water waste with exclusive Elix® technology.
- Experience More www.millipore.com/ultrapure



PEAKS of Interest

Daniel W. Armstrong to Receive ACS Award in Separation Science and Technology 2014

Daniel W. Armstrong, a professor of chemistry and biochemistry at the University of Texas Arlington in Arlington, Texas, is to receive the 2014 American Chemical Society Award in Separations Science and Technology. This is the second time this year that the ACS has recognized the work of Armstrong — in August last year he was selected to join the ACS Class of Fellows 2013.

The award, sponsored by Waters Corporation, recognizes “the development of novel applications with major impacts and/or the practical implementations of modern advancements in the field of separation science and technology.” It will be presented to Armstrong on March 17 at the ACS national conference in Dallas, Texas.

Armstrong is thought of as a leader in his field and is considered to be the “father” of pseudo-phase separations — a type of chromatography that lowers costs, volatility, and toxicity while providing higher selectivity than other analytical methods. His commercialized innovations are applied across the drug development, environmental analysis, and petrochemical industries.

In his career, Armstrong has developed more than 30 different chromatography columns — a gas chromatography column that he developed is now part of the Rosetta mission of the European Space Agency exploring the composition of comets in space.

Purdue Research Foundation Invests in Two Start-up Companies

The Purdue Research Foundation (West Lafayette, Indiana) has awarded \$150,000 in investment to two companies that have licensed Purdue University Innovations: Telos Discovery Systems Inc., and bioVidria Inc. The investment is from the Emerging Innovations Fund, which was set up in 2008 to accelerate the commercialization of early stage technologies in the Purdue community.

The company bioVidria Inc., received \$75,000. It is based on the work of Mary J. Wirth, the W. Brooks Fortune Distinguished Professor of Analytical Chemistry in the Department of Chemistry. The company focuses on developing new materials for chromatography, specifically protein analysis, for the pharmaceutical and agricultural industries.

Telos Discovery Systems Inc., also received \$75,000. The company is a boutique equipment manufacturer that supports basic and discovery medical research. The company is based on the work of Joseph Garner, a former Purdue associate professor in the Department of Animal Sciences.

Elizabeth Hart-Wells, assistant vice president and director of the Purdue Office of Technology Commercialization said: “The fund makes seed investments and early-stage investments in the form of convertible debt and/or warrants, with a goal of accelerating the commercialization of early-stage innovations.” ■

LC|GCtv New videos from LCGC



KOEN SANDRA ON ANALYZING BIOPHARMACEUTICALS WITH LC-MS

By 2020, more than 50% of approved drugs are expected to be protein biopharmaceuticals. Koen Sandra from the Research Institute for Chromatography in Belgium discusses the challenges of characterizing protein biopharmaceuticals and why LC-MS is now indispensable for characterizing this class of drugs.

Other recent LCGC TV interviews include:

- Dwight Stoll of Gustavus Adolphus College on 2D LC
- Mary Ellen McNally of DuPont Crop Protection on the proper use of an internal standard
- Michal Holcapek on trends in LC-MS
- Frederic Lynen on analyzing drugs and metabolites in biological matrices

Visit <http://www.learnpharmascience.com/lcgc> to see these videos and more.

IS THERE AN APP FOR THAT?

HIGHLIGHTS OF APPS FOR SEPARATION SCIENCE

4G 100% 11:59 AM

APP NAME: HPLC Troubleshooting Guide

OFFERED BY: Thermo Scientific

PLATFORMS: Android, iPhone (optimized for iPhone 5), iPad, iPod Touch

WHAT IT DOES: This app diagnoses various potential issues with HPLC analysis. Solutions range from general nature issues that are common to standard instrument setups from various manufacturers to troubleshooting special technologies such as nano LC and fluorescence, charged aerosol, and electrochemical detection. The app also offers interactive functions to connect to Thermo Scientific social media, recommend the app to a friend, and contact the company for further support. By tapping the Share function, troubleshooting solutions can be shared among teams.

COST: Free





Now for my next trick: Essential Macromolecular Characterization™ *without* SEC-MALS

If you're not using Wyatt Technology's Multi-Angle Light Scattering detectors coupled with Size Exclusion Chromatography (SEC-MALS), Field Flow Fractionation (FFF-MALS), or Composition Gradient (CG-MALS) instruments to make your essential macromolecular characterization measurements, or the DynaPro® Dynamic Light Scattering (DLS) Plate Reader for high-throughput screening, you must really believe in magic! When you add a Wyatt MALS, DLS or mobility detector to your laboratory, you'll be able to determine *absolute* molar masses and sizes, quantify biomolecular interactions and charge, and run automated screens of size, aggregation and stability, completely independent of the typical legerdemain of ancient analytical techniques. Your colleagues will definitely be amazed.



DAWN® HELEOS® Multi-angle light scattering for characterizing macromolecular mass, size and conformation. No other MALS detector comes close in range or sensitivity.



Optilab® T-rEX™ Differential refractometer for on-line concentrations. Unparalleled measurement and temperature ranges.



Mobius® Electrophoretic mobility from proteins to nanoparticles. Uniquely sensitive and versatile for all types of solvents.




Eclipse™ Column-free, field-flow fractionation of macromolecules and nanoparticles in solution. One system covers all your separation needs.



DynaPro® Plate Reader Automated dynamic light scattering for proteins and nanoparticles, *in-situ* in standard micro-well plates. Now with a camera!



WYATT
TECHNOLOGY
www.wyatt.com



COLUMN WATCH

In-Capillary Sample Concentration in CE

In this installment, guest authors and capillary electrophoresis (CE) experts Breadmore and Sanger-van de Griend bring a practical perspective to enhancing the sensitivity of CE. They compare stacking and sweeping injection methods and provide practical suggestions about how to enhance sensitivity for charged and neutral compounds.

Michael C. Breadmore and Cari E. Sanger-van de Griend are the guest authors this month. **Ronald E. Majors** is the editor of *Column Watch*.

One of the most frequently cited concerns in capillary electrophoresis (CE) is that the concentration limits of detection are inferior to what can be achieved with liquid chromatography (LC) (1). Extensive research has been carried out to overcome this limitation and now there are many different ways to increase sensitivity by controlling the way the sample is injected, how it is prepared, and how the targets are separated. After all, the concentration in the sample vial does not need to be the same as the concentration needed for detection if we can concentrate on-line. Early researchers in the field diluted samples in the electrolyte (parallel to dissolving or diluting in eluent in LC); however, doing so dilutes the target analytes, further reducing sensitivity. A better approach is to use the sample straight, which means that the composition of the sample plug is different from the composition of the electrolyte. Although intuitively this may not seem like a smart idea, depending on the exact analyte and conditions there may be a significant difference in their migration behavior in the sample zone compared to their migration behavior in the electrolyte. Very early researchers realized that this approach could enhance sensitivity. A brief search of the literature yields many excellent research papers and reviews (1–5), with different kinds of concentration mechanisms occurring because of different analyte, matrix, and electrolyte properties. Table I gives an overview of the most frequently used concentration techniques; as can be seen there are quite a number of them. It is worthwhile noting that these are the most common single concentration approaches used. It is possible to

combine multiple mechanisms to produce even more powerful approaches; however, these are more complex. Most reviews discuss the various techniques from the mechanistic point of view. Here, we'd like to consider these options from a practical point of view and ask the question: I have XYZ in my sample; how do I improve my sensitivity? After all, in real life one usually cannot pick all conditions and sample matrices to suit a theoretically useful technique; instead, one has to work from the situation at hand.

Velocity Differences

First of all, the fundamental premise of concentration in CE comes back to creating a different velocity of the analytes of interest in the sample zone compared to their velocity in the background electrolyte (BGE) (6). There are many different techniques based on a change in velocity. Each technique has its own abbreviation, but the only difference may be a change in direction of the electroosmotic flow (EOF). These techniques fall into two groups: stacking techniques and sweeping techniques. If the velocity in the sample zone is higher than in the BGE (Figure 1b), the technique is generally called *stacking* (7). When stacking, the analytes quickly migrate through the sample zone to the electrolyte boundary where they slow down and "stack" together on the boundary with an increase in concentration and, thus, sensitivity (Figure 1b). In so-called *sweeping* techniques, in contrast, the velocity in the sample zone is slower than the velocity in the electrolyte. As the boundary between the electrolyte and sample moves through the sample matrix, analytes are concentrated around this

JUST GOOD
ENOUGH...
ISN'T FOR US

Innovating GC
solutions for
challenging applications.

ENGINEERING ADVANTAGE

**YOU CAN SEE
OUR THINKING**

*Downtime from rotors?
Come see us at
PITTCON Booth 639*



See more of our thinking and the Advantages it delivers.
Call +1.303.794.2611 or email enquiry@usa.norgren.com
<http://programs.norgren.com/NorgrenAFP>

Table I: In-capillary sample concentration in CE

Stacking	The analyte velocity in the sample zone is faster than the analyte velocity in the BGE
Electric Field Strength Induced Velocity Differences	
Field-amplified sample stacking (FASS)	Sample has low conductivity (<1/10 BGE), so the local field strength is high and local velocities are high. The analytes slow down and are stacked when entering the BGE.
Field-amplified sample injection (FASI)	Sample in a low-conductivity matrix is injected electrokinetically. The local field strength is high, increasing the amount injected. If the EOF can be suppressed, long injection times are feasible. Injection is discriminative.
FASI with solvent plug	Same as FASI, but a solvent plug is injected with the sample for low conductivity and high local field strength.
Large volume sample stacking (LVSS)	A large sample plug is injected. The sample matrix is continuously removed by applying a reversed voltage during stacking.
Large volume sample stacking with an EOF pump (LVSEP)	A large sample plug is injected. The sample matrix is continuously removed by the reversed or suppressed EOF in the sample zone (compared to the EOF in the BGE).
Micellar electrokinetic chromatography (MEKC) stacking	Low concentration of pseudo-stationary phase (PSP) moves rapidly through low conductivity sample (<1/10 BGE). PSP with analyte stacks at boundary between sample and electrolyte.
pH-Induced Velocity Differences	
Dynamic pH-junction	The analytes have a higher mobility in the sample zone and will be stacked at the boundary with the BGE
pH-mediated FASS	The sample is injected between a strong base and strong acid. When the voltage is applied, the OH ⁻ and H ⁺ migrate through the sample zone and neutralize, thus reducing the conductivity. Stacking can then take place.
Sweeping	
Sweeping	The analyte velocity in the sample zone is slower than the analyte velocity in the BGE. Sample without PSP is injected onto a system with a PSP BGE. The analytes are concentrated at the front of the PSP zone moving through the sample zone. Conductivity of sample equivalent to BGE.
Micelle to solvent stacking (MSS) and analyte focusing by micelle collapse (AFMC)	The sample is mixed with micelles, the BGE is micelle-free. The micelles collapse when meeting the solvent or the BGE. The analytes are concentrated at the boundary.
Isotachopheresis	
Transient isotachopheresis (tITP)	The sample is injected between a higher mobility co-ion (leading electrolyte) and a lower mobility terminating electrolyte. When the voltage is applied, the analytes sort in order of mobilities. The concentrations in the analyte zones depend on the concentration of the LE (Kohlrausch regulation function).
In-Capillary SPE	
In-capillary solid-phase extraction (SPE)	A small SPE column is manufactured inside the capillary. Analyte is extracted through conventional chromatographic principles and eluted in solvent for electrophoretic analysis. Easily automated, attractive for low-volume samples.
Liquid-Phase Microextraction	
Single-drop microextraction (SDME)	Analytes are extracted into a small organic drop suspended at the tip of the capillary. Contents of the drop are then analyzed by electrophoresis. Can be performed in commercial instrumentation with no modifications.
In-vial supported liquid membrane (SLM)	Organic liquid is held in a thin membrane. Analytes are extracted from sample through organic phase into aqueous acceptor. Can be performed in-vial.
In-vial electrokinetic microextraction (EME)	Similar to in-vial SLM, except a voltage is applied across the membrane to enhance transport of charged species through the organic phase.

moving boundary and are considered to be “swept” together with an increase in concentration and, thus, sensitivity (Figure 1c) (8).

The Easiest Way: FASS

So if we want to create a concentrating effect inside the capillary, we have to ensure that the analyte has a difference in velocity between the sample zone and BGE. In an electric field, the velocity of an analyte is directly proportional to its electrophoretic mobility, which is governed by charge, size, and the electric field strength. This means that we need to know the analyte’s physical-chemical properties. The first

question in the decision process (see Figure 2) is as follows: Is the analyte charged? (And yes, whether the species is anionic or cationic is important, but only for determining the conditions to use when implementing the method; this does not matter when choosing the ideal approach for concentration.) If the analyte is neutral, then you need to use affinity approaches to achieve concentration, which are discussed later on.

The second question is as follows: Does the sample matrix have a lower conductivity than the BGE? If yes, then *field-amplified sample stacking* (FASS) may well work (Figure 1b). This is the simplest approach because the lower

conductivity of the sample zone automatically means that the field strength in the sample zone is higher than in the BGE. If the mobility of the analyte remains the same, then stacking will occur as a higher field strength results in a higher velocity without any other active manipulation (there is nothing you can do about it except to remove the conductivity difference). This is the reason that it is preferred to have the sample dissolved or diluted in something that has a lower conductivity than the BGE, including diluted BGE, pure water, or organic solvents. Ethanol, dimethyl sulfoxide (DMSO), and acetone can give even better results

GPC Innovations from the Experts

Turn up the Heat on your GPC Analysis!


With the new EcoSEC® High Temperature GPC System and TSKgel® Columns

The EcoSEC High Temperature GPC System incorporates the proven design and technology used in our renowned ambient system with new TSKgel H series columns, designed for durability up to 220 °C.

This combination delivers reliable, reproducible results: up to the challenge of your high temperature analyses.

Learn more by visiting the Tosoh booth #807 at Spring ACS and visit our website to register for the Advanced Polymer Characterization Workshop.

*Tuesday, March 18 from 12:30 PM – 3:00 PM
Expo Halls C/D, Workshop Room 2.*



TOSOH BIOSCIENCE LLC • Customer service: 866-527-3587 • Technical service: 800-366-4875, option #4

Tosoh Bioscience, TSKgel and EcoSEC are registered trademarks of Tosoh Corporation.

TOSOH BIOSCIENCE
www.tosohbioscience.com






Figure 1: Concentrating in CE. Creating velocity differences of the analyte in the sample zone versus the velocity of the analyte in the BGE. (a) Initial situation after injection and before applying the voltage. (b) Stacking, the analyte has a higher velocity in the sample zone than in the BGE. (c) Sweeping, the analyte has a lower velocity in the sample zone than in the BGE.

How Much Is Injected?

How much sample you should inject in CE depends on your injection method.

In hydrodynamic injection, the injection volume is determined by the following equation:

$$V = \frac{\Delta p d^4 \pi t}{128 \eta L}$$

The volume (V) of sample injected depends on the injection time (t) and pressure difference (Δp), the viscosity (η), the capillary diameter (d), and total capillary length (L).

In electrokinetic injection, the injection volume is determined by the following equation:

$$Q = \frac{(\mu_e + \mu_{eo}) E \pi d^2 C t}{4}$$

The amount (Q) injected of a certain analyte depends on the analyte's effective mobility ($\mu_e + \mu_{eo}$) and concentration (C), the capillary diameter (d), the applied electric field (E), and the injection time (t). Electrokinetic injection is a selective or discriminative injection mode as analytes with different mobilities are injected in different amounts.

than pure water (9). As a consequence of the stacking effect, a larger volume than the usually mentioned 1% of the capillary volume can be injected. As a rule of thumb, injection volumes up to 5% of the total capillary volume are feasible with FASS, with enhancements in sensitivity typically less than about 50 achieved. Larger volumes do not produce further improvements because of mismatches in EOF, which causes hydrodynamic dispersion, a loss of efficiency (10,11), and decreased reliability because of the occurrence of bubbles.

When a Higher Concentration Is Needed

The sidebar "How Much Is Injected" provides some help in calculating the injected volumes or amounts. Sometimes a greater increase in concentration is required than can be achieved in FASS. This is easiest if the sample matrix has a controlled and constant composition. In such a case, *field-amplified sample injection* (FASI) (12) (also called *field-enhanced sample injection* [FESI]) and *large-volume sample stacking* (LVSS) can be easily used. The only difference between FASS and FASI is the injection — for FASI, electrokinetic injection is used. Because the amount injected is strongly dependent on the local field strength, the sample matrix conductivity should be constant. A variable sample matrix composition will result in variable amounts injected, which is why for all practical purposes FASI requires a controlled sample composition. The same is true for LVSS for similar reasons. In LVSS, the sample matrix is removed either through polarity switching (denoted *LVSS*) (10) or EOF pumping (denoted *large-volume sample stacking with an EOF pump* [LVSEPI]) (13).

If the matrix varies, local field strengths vary and timing becomes an issue, resulting in variable concentration effects and migration times. To some extent these problems can be compensated for with an internal standard, but method performance (detection limit, repeatability, and so on) will be compromised so it is better to address the cause of the problem, rather than doing a patch-up job with duct tape and silly putty.



SELECTRA[®] C18 HPLC Columns

- Good choice for easily converting current C18 methods
- Suitable for traditional reverse phase analyses
- High purity, type B spherical silica
- Excellent Carbon loading
- Fully endcapped
- Conforms to USP L1



NEW!



800.385.3153 www.unitedchem.com




Figure 2: In-capillary sample concentration decision tree for charged analytes.




Figure 3: Isotachopheresis (ITP), (a) initial conditions with injection of sample between the leading and terminating electrolytes, (b) isotachopheretic zones at steady state.

If the matrix composition is variable, we first look again at the analyte: Is it a weak base or acid? If yes, stacking with a *dynamic pH junction* (6,14), also

called a *moving neutralization boundary* (15), is feasible. Here, the BGE and sample composition are selected such that there is a difference in charge,

because a charge difference will also result in different velocities between the sample and BGE zone. It is important to note that it is not necessary to completely reverse the charge; a simple change from charged to neutral may be sufficient if the system is designed correctly. If the answer is no, sweeping techniques may be appropriate. In sweeping, analyte affinity for a pseudo-stationary phase (PSP) is key. When the PSP migrates through the sample zone, analyte molecules are concentrated at the front of the PSP as it moves through the sample zone (Figure 1c) (8). If the PSP is a part of the BGE, such as in micellar electrokinetic chromatography (MEKC), it is called *sweeping*. If the PSP is in the sample and the BGE is PSP-free, then when the micelles migrate into a PSP-free zone, the concentration of the surfactant reduces to below the critical micelle concentration and the micelles collapse, releasing the analytes. This was first introduced for neutral species as *analyte focusing by micelle collapse* (AFMC), but has been applied to charged species with organic solvents to collapse the micelles and is then called *micelle to solvent stacking* (MSS). It is important to note that the mechanism of sweeping and MSS still function with a variable matrix, but there may be some resultant variability in the final separation. Nevertheless, similar to isotachopheresis (ITP), discussed below, these approaches are very useful for samples with a variable matrix.

High-Conductivity Sample Matrix

If the sample matrix has a high conductivity, it is important to consider off-line cleanup. From the discussion above and from the decision tree (Figure 2), it can be seen that with a lower conducting sample matrix, many approaches can be used for successful sensitivity enhancement. If the sample matrix has a high level of conductivity but off-line cleanup is not feasible, there are some options available. These are based on creating velocity differences by sweeping or MSS, by manipulating the charge by pH-mediated techniques (such as pH-mediated FASS or dynamic pH junction), or by ITP.


Of these, *pH-mediated FASS* is an interesting and unique approach that induces a FASS effect in an indirect manner (16). In this approach, sample is injected electrokinetically, and for anions, the electrolyte would contain a weak base, like ammonium, such that during injection the anions migrate into a zone where the counterion is ammonium (not sodium as it may have been in the sample). After injection of the sample, there is a short electrokinetic injection of base (or acid for cations). The base migrates into the sample zone (hydroxide and hydronium have a high mobility) and titrates the buffer counterion into water, thereby producing neutral species and reducing the conductivity of the sample zone. With this change in conductivity, there is now a higher electric field over the sample and FASS can occur. Although this is a very clever approach, it is not easy to optimize the conditions to achieve the outcome required.

Transient ITP

Perhaps the most attractive approach for charged species in samples with a high conductivity matrix is to use *isotachopheresis* (ITP) (Figure 3). ITP is one of the oldest forms of CE and probably the most intimidating — and also poorly understood by many experienced practitioners of CE. In ITP, the sample is sandwiched between a leading electrolyte and a terminating electrolyte. The leading electrolyte has an ion with the same charge as the analyte but with a higher mobility. The terminating electrolyte has an ion with the same charge as the analyte but a lower mobility. The sample analytes sort in migration order and are concentrated to obtain the same conductivity in each analyte zone. The mobility and concentration of the leading electrolyte determine the concentrating effect on the analytes, according to the Kohlrausch regulation function, and the mobility window between the leading and terminating ions provides a zone in which the analyte range can be selectively tuned. The unique mechanism of ITP makes it the great equalizer. Macrocomponents are diluted and trace components are concentrated, making ITP exceptionally useful as a stacking technique.

The simplest way to use ITP is to perform transient ITP (tITP) (17). Here, a short ITP stage is induced before the CE separation. Transient ITP can be done by injecting small volumes of leading and terminating electrolytes on either side of the sample; alternatively, either the leading or termination co-ion comes from the BGE (not both) and the other comes from the sample or from an additionally injected plug. Transient ITP is particularly suited to salty samples — especially those that

contain lots of sodium chloride as the chloride is a very high mobility ion that can function as a leader for almost all anions. Sodium is a mid-range mobility cation so it cannot be used for all cations, but it typically has a higher mobility than many organic cations. It is also worth briefly mentioning *pseudo-ITP* (pITP), in which the terminating zone is created not from a low mobility ion, but from a solvent added to the sample in which there is a lower conductivity (18). This approach arose out of work



Polymicro Capillary Tubing

TRUST THE INDUSTRY LEADER > TO BRING THE BEST SOLUTIONS TOGETHER

Polymicro Technologies, a subsidiary of Molex, offers expertise and true innovation in capillary tubing — including our new 1/32" OD tubing for LC applications. From custom sizes and geometries, to laser machining and precision cleaving, Polymicro capillary tubing ensures your critical success. And when custom assemblies are required, you can be confident that you've got the industry leader at work for you, delivering only the highest quality solutions for your application. Count on us for all of your capillary tubing needs.

See us at Pittcon, Booth 2224

Polymicro
TECHNOLOGIES

A Subsidiary of **molex**

www.molex.com/lcgc/polymicro




Figure 4: In-capillary sample concentration decision tree for neutral analytes.

where the supernatant of biological samples after deproteination with acetonitrile was injected into the CE, and it is therefore a very useful technique for plasma and serum samples because they require very simple sample processing.

The Last Resort: In-Capillary and In-Vial Extraction

If stacking and sweeping approaches do not work, then the only real alternative is to consider the use of an on-line or in-capillary extraction. Such extractions have traditionally been performed by making a small chromatographic bed (19) at the tip of the capillary (most easily done with larger-internal-diameter tubing and particles that have a larger diameter than the internal diameter of the separation capillary) and allowing analytes to be extracted chromatographically before electrophoretic separation. The technique resembles on-line *solid-phase extraction* (SPE) and includes possibilities for immunoaffinity extraction as well (20). While this is a potentially very attractive approach, there are a number of limitations including analysis time, repeatability, and the lack

of commercially available capillaries. Therefore, it is really only useful for researchers.

The most common off-line alternative to SPE is *liquid-liquid extraction* (LLE) (21), which has only very recently been reported coupled directly with CE through developments in *single-drop microextraction* (SDME). Small individual drops of organic and organic-aqueous drops can be created on the tip of the capillary in a commercial CE unit (22). This allows for easy automation, with very good analytical performance; however, this is not a straightforward approach at this point in time because a slight back pressure is required to stabilize the drop. A more attractive approach is the implementation of the technique in a vial. This can be done with a water-immiscible ionic liquid (23), and excellent results have been obtained from urine and serum; this is an area very much still in its infancy, and it is unclear yet whether this technique will be attractive. The most promising approach to appear recently is the development of in-vial *supported liquid membranes* (SLM) (24) and *elec-*

tromembrane extraction (EME) (25). Still very much in the research phase, these techniques involve creating a custom vial in which the sample is placed below a membrane impregnated with a suitable extraction solvent with an aqueous acceptor phase placed on top of the membrane. When placed in a commercial CE system, the capillary samples directly from the surface of the membrane, which considerably decreases the time required to reach localized equilibrium (26). This approach has been demonstrated for inorganic ions and a few drugs from various biological fluids. If it can be transitioned into a commercial vial, it will significantly aid the ability of CE to easily handle complex and high conductivity samples.

How to Deal with Neutral Species

The options for uncharged analytes are fewer than those for charged analytes. For uncharged analytes, all in-capillary approaches are based on some kind of affinity interaction (27). This is because neutral species are uncharged and do not have their own electrophoretic mobility and, thus, cannot be concentrated by approaches developed for charged compounds. In contrast, many of the approaches developed for neutral analytes can also be applied to charged analytes and some were already discussed in detail above. Figure 4 shows the decision tree for uncharged analytes.

The first approaches to concentrate neutral species were based on stacking: By using a low conductivity sample and affinity between the analyte and a PSP, it is possible to rapidly transport analyte molecules from the sample zone to the sample or electrolyte boundary where they concentrate. There is an abundance of slightly different approaches here depending on the direction of the micelle and the polarity, but they all rely on stacking (not sweeping) for concentration (27).

Today, however, the method of choice for concentrating neutral analytes is sweeping as this can provide greater enhancements than stacking and is applicable to samples with a low, high, or variable sample matrix (8,28). The challenge with sweeping is that for good enhancement a high affinity with the PSP is required (the zone length

after sweeping is related to the affinity interaction), but this approach is bad for separation as all of the analytes will migrate with the PSP. Therefore, sweeping requires careful optimization to ensure that sufficient sensitivity enhancements can be obtained. Sweeping is also attractive because it is applicable to high conductivity samples — the mechanism of concentration does not depend on the conductivity of the sample. In fact, when sweeping is combined with electrokinetic injection, multiple capillary volumes of sample can be injected by the EOF with the analytes concentrated and held at the front of the micelle zone (29). This approach can provide very significant enhancements, but the trade-off is time — typically, very long injections (approaching an hour) are required.

The Chicken or the Egg — Which Part of Your Method to Consider First

The typical approach used by most workers in the field is to develop separation conditions to resolve the target analytes from matrix components, and to then consider whether or not any on-line concentration techniques can be used. The justification here is twofold; it doesn't really matter what sensitivity you have if you can't separate the analytes (with increased prevalence of mass spectrometers this may become arguable) and you can always resort to off-line treatment of your sample to improve the sensitivity. This latter point in particular is worth some consideration because it is not easy, for example, to concentrate by more than 100 — this is the equivalent of taking 10 mL of sample, and extracting the analytes into a volume of 100 μ L. A volume of less than 100 μ L can be difficult to handle and introduces other errors. Stacking and sweeping approaches are getting closer to being able to electrokinetically extract from 100 μ L into 1–5 nL (100,000-fold), but this is only possible for certain sample types and is only possible if the separation method is compatible with these approaches.

You don't need a chicken or an egg — you need both. Therefore, it would be prudent to consider the requirements of sensitivity during method development and ensure that the developed method is compatible with at least some of the concentration approaches. For example, if


you know that you will need nanogram-per-liter detection limits or lower, then that really does place restrictions on the concentration approaches that can be used and implicates which methods you need to use. This, in turn, defines the electrolyte compositions that may be more ideal and compatible.

Other Practical Limitations: pH, pI, Stability, and Analyte Size

In the discussion of sensitivity enhancement so far, we have not indicated

whether any technique is typically more or less suitable for a certain class of molecules, other than via physicochemical properties like weak acid and bases, charged or uncharged analytes, and high or low affinity for a pseudo-stationary phase. This is deliberate because it is the physicochemical properties that govern which method is suitable. A dynamic pH junction (essentially a step pH gradient) is useful for proteins with an isoelectric point (pI) between the pH of the sample and the electrolyte, not because of size,

RegisPack® • RegisCell® • Wheelk-O®1 • RegisPack CLA-1™



Excellent Selectivity

High Quality

Superior Separations

SFC Separation Services

Free Chiral Screening


Add power to your chiral separations

Regis Technologies' complete line of chiral columns separates over 90% of all chiral compounds. The broad selection of CSPs includes original Pirkle-Type columns, coated polysaccharides and crown ether CSPs. Our FREE chiral screening determines the best chiral column for your enantiomer analysis. We are the chiral chromatography experts.

It all adds up to the Power of Regis.

847.967.6000

www.RegisTech.com



CHROMATOGRAPHY GMP MANUFACTURING SFC SERVICES PHARMA SERVICES API FINE CHEMICALS

REGIS TECHNOLOGIES CHIRAL CHROMATOGRAPHY




Figure 5: Sweeping and FASS plus sweeping of standard mixture of hypolipidemic drugs. Capillary: 78 cm \times 50 μ m LPA-coated capillary; BGE: 20 mM ammonium bicarbonate (pH 8.50) containing 50 mM SDS; voltage: -28 kV. Sweeping: hydrodynamic injection of 10 μ g/mL drug mixture in ammonium bicarbonate with conductivity similar to BGE at 50 mbar for 40 s. FASS plus sweeping: hydrodynamic injection of 1 μ g/mL drug mixture in Milli-Q water at 50 mbar for 80 s. Peaks: A = atorvastatin, F = fluvastatin, G = gemfibrozil, L = lovastatin, M = mevastatin, P = pravastatin, R = rosuvastatin, S = simvastatin. Adapted from reference 30 with permission from Elsevier.

Good CE Injection Practice

Injection times should be

- As short as possible
- But at least 3 s to compensate for injection pressure or voltage variability

Injection plug length less than 1–2% of the total capillary length

Constant, precise temperature for capillary and vials

Constant level at outlet vial during injection, no waste vial

Sample vials same levels

Background electrolyte (BGE) inlet and outlet vials same levels

Dip capillary inlet in clean water or BGE vial after injection

Inject BGE or water plug after sample injection

Burn off polyimide capillary ends (not for certain coated capillaries)

Cut the capillary ends straight

Sample matrix of samples and standards should match

When using an internal standard:

- Reduce injection volume variability
- Reduce injected amount variability
- Correct for migration time variability

When using a stacking technique:

- Dilute or dissolve sample in water or diluted BGE
- Use literature for inspiration to gain sensitivity

When using electrokinetic injection:

- Selective
- Only one injection per vial
- Variable ionic strengths influence injected amount
- Amount injected related to d^2 and L

When using hydrodynamic injection:

- Generally more precise than electrokinetic injection
- Not appropriate for viscous applications
- Multiple injections per vial feasible
- Volume injected related to d^4 and L

but because of their pI values. Of course there are other practical limitations that need to be considered for all molecules, but have some greater relevance for larger molecules. Knowledge about other parameters such as stability or adsorption should always be included in the decision process. Limiting for all concentration efforts is the solubility of the analyte — low solubility situations should be avoided (such as proteins at their pI) unless a solubility-enhancing solvent or compound is added. The sidebar “Good CE Injection Practice” lists some other practical tips to improve method performance.

What’s Good for the Goose, Isn’t Always Good for the Gander

As can be understood from the discussions above, opportunities for sensitivity enhancement through injection are strongly influenced by other method parameters, such as the BGE composition. A higher BGE concentration, for example, makes it easier to increase the conductivity difference with the sample, but this causes higher currents, which may lead to greater method instability. Although very high sensitivity enhancements have been published with stacking techniques, in the end it always comes down to compromising between robustness and resolution (that is, separation plus efficiency) versus sensitivity. Depending on the required use one has to make choices and take the consequences.

A Scot designed the
much loved Buick.
We also design the
best gas generators.



David Dunbar Buick, was an inventive Scot who now has his name stamped on over 25 million cars in the USA. But Scottish inventiveness didn't stop there...

Today, Peak Scientific gas generators also proudly represent Scottish innovation. Manufactured to ever higher standards and with new features appearing regularly, Peak's gas generator range offers innovation, reliability, and the assurance of peak performance.

Add to this, Peak [**Protected**] world class after sales care and you can see why Peak generators are parked in some of the most sophisticated labs in the world.

FIND OUT MORE

Call: +1 866 647 1649

Email: marketing@peakscientific.com

www.peakscientific.com

PEAK 
SCIENTIFIC

Expect Innovation

Stacking or Sweeping — Or Both?

The right in-capillary sample concentration method to use with CE depends on the properties of the analytes. With a target set of analytes with diverse properties, finding a suitable method can be complex. An example of this is shown in Figure 5. This is a mixture of eight hypolipidemics with a diverse range of properties including mobility, pK_a , and lipophilicity. This latter property is a problem because pravastatin, gemfibrozil, and rosuvastatin show poor interaction with the PSP and thus are not swept well. The solution in this case was to perform sweeping under FASS conditions — sweeping the high lipophilic analytes and stacking the low lipophilic ones. The enhancements are modest (40–80-fold) but this nevertheless shows the potential of using multiple mechanisms simultaneously to improve analyte coverage. This is different than using sequential approaches such as dynamic pH junction—sweeping, which uses one method to concentrate and then a second method to sharpen the peaks further and improve efficiency. These combined approaches are more complex and beyond the scope of the article, but there really is no limit to the combinations of methods that can be used.

The Final Word

The right approach to sample concentration for CE is one that will meet your needs. Determining this isn't a simple matter, and in many cases there may be more than one possible approach that would be suitable. Literature can be useful to identify what may be more appropriate, but ultimately it will come down to experimenting until the appropriate conditions are found. The best advice is to start with a simple approach and only add complexity if it is absolutely needed.

Acknowledgment

The sidebar "Good CE Injection Practice" was reproduced with permission, Cari Sanger, CESolutions 5. Injection precision and sensitivity, Separation Science, <http://www.sepscience.com/Techniques/CE/1352-/CE-Solutions-5-Injection-Precision-and-Sensitivity>.

References

- M.C. Breadmore, *J. Chromatogr. A* **1221**, 42–55 (2012).
- A. Slampova, Z. Mala, P. Pantuckova, P. Gebauer, and P. Bocek, *Electrophoresis* **34**, 3–18 (2013).
- M.C. Breadmore, A.I. Shallan, H.R. Rabanes, D. Gstoettenmayr, A. Syazwani Abdul Keyon, A. Gaspar, M. Dawod, and J.P. Quirino, *Electrophoresis* **34**, 29–54 (2013).
- F. Kitagawa and K. Otsuka, *J. Chromatogr. A* (2013), <http://dx.doi.org/10.1016/j.chroma.2013.10.066>.
- Y. Wen, J. Li, J. Ma, and L. Chen, *Electrophoresis* **33**, 2933–2952 (2012).
- P. Britz-McKibbin, G.M. Bebault, and D.D.Y. Chen, *Anal. Chem.* **72**, 1729–1735 (2000).
- R.-L. Chien and D.S. Burgi, *Anal. Chem.* **64**, 489A–496A (1992).
- J.P. Quirino and S. Terabe, *Anal. Chem.* **71**, 1638–1644 (1999).
- Z.K. Shihabi, *J. Chromatogr. A* **1066**, 205–210 (2005).
- R.-L. Chien and D.S. Burgi, *Anal. Chem.* **64**, 1046–1050 (1992).
- C. Huhn and U. Pyell, *J. Chromatogr. A* **1217**, 4476–4486 (2010).
- R.-L. Chien and D.S. Burgi, *J. Chromatogr.* **559**, 141–152 (1991).
- Y. He and H.K. Lee, *Anal. Chem.* **71**, 995–1001 (1999).
- P. Britz-McKibbin and D.D.Y. Chen, *Anal. Chem.* **72**, 1242–1252 (2000).
- C.-X. Cao, L.-Y. Fan, and W. Zhang, *Analyst* **133**, 1139–1157 (2008).
- S. Park and C.E. Lunte, *J. Microcol. Sep.* **10**, 511–517 (1998).
- F. Foret, E. Szoko, and B.L. Karger, *J. Chromatogr.* **608**, 3–12 (1992).
- Z.K. Shihabi, *J. Chromatogr.* **652**, 471–475 (1993).
- R. Ramautar, G.W. Somsen, and G.J. de Jong, *Electrophoresis* **35**, 128–137 (2014).
- N.A. Guzman, T. Blanc, and T.M. Phillips, *Electrophoresis* **29**, 3259–3278 (2008).
- S. Pedersen-Bjergaard and K.E. Rasmussen, *Anal. Chem.* **71**, 2650–2656 (1999).
- K. Choi, J.K. Su, G.J. Yoo, O.J. Yong, J.S. Kim, and S.C. Doo, *Anal. Chem.* **81**, 225–230 (2009).
- M.C. Breadmore, *J. Chromatogr. A* **1218**, 1347–1352 (2011).
- S. Palmarsdottir, E. Thordarson, L.E. Edholm, J.. Jonsson, and L. Mathiasson, *Anal. Chem.* **69**, 1732–1737 (1997).
- S. Pedersen-Bjergaard and K.E. Rasmussen, *J. Chromatogr. A* **1109**, 183–190 (2006).
- P. Pantuekova, P. Kuban, and P. Boek, *J. Chromatogr. A* **1299**, 33–39 (2013).
- A.T. Aranas, A.M. Guidote, Jr., and J.P. Quirino, *Anal. Bioanal. Chem.* **394**, 175–185 (2009).
- J.P. Quirino and S. Terabe, *Science* **282**, 465–468 (1998).
- J. Palmer, N.J. Munro, and J.P. Landers, *Anal. Chem.* **71**, 1679–1687 (1999).
- M. Dawod, M.C. Breadmore et al., *J. Chromatogr. A* **1217**, 386–393 (2010).

This month's guest coauthors:

Michael C.

Breadmore is an ARC Future Fellow (2013–2017) at the Australian Centre for Research on Separation Science (ACROSS) based at the University of Tasmania. His research interests lie in the application of capillary and microchip electrophoresis to the trace analysis of environmental, clinical, and forensic samples, as well as in methods for the low-cost fabrication of microfluidic devices.



Cari E. Sanger-van

de Griend is a consultant and trainer at Kantisto BV, a pharmaceutical analytical consultancy she founded after 20 years in industry. She is a globally recognized expert on the implementation of CE techniques and a specialist on pharmaceutical analysis. She is also an Associate Professor in Analytical Pharmaceutical Chemistry at Uppsala University in Sweden and an Adjunct Senior Lecturer at the Australian Centre for Research on Separation Science (ACROSS), University of Tasmania.



The editor of Column Watch:

Ronald E. Majors

"Column Watch" Editor Ronald E. Majors is an analytical consultant and is a member of LCGC's editorial advisory board. Direct correspondence about this column to lcgcedit@lcgcmag.com.



For more information on this topic,
please visit
www.chromatographyonline.com/majors

John McFee's high tech hospital clinic

Dried blood spot (DBS) sampling is an emerging technology in the clinical and pharmaceutical laboratory, offering easy sample collection, transport and storage. Our revolutionary DBS Autosampler maintains the integrity of the sample through automation, offering time and cost savings. Innovative patented Flow-through desorption technology (FTD™)* eliminates tedious punching and costly robotics. Automation of the entire workflow for DBS analysis in minutes, providing maximum sensitivity without any manual intervention.


We invite existing and new partners in the clinical and pharmaceutical areas to an exclusive preview at our booth at MSACL (#32), Pittcon (#1839) and analytica (#A1.407).

- **Automated workflow – rapid results**
- **Maximum sensitivity**
- **Minimal sample transport and storage costs**
- **Ease of sample collection**

Spark Holland B.V.
P.O. box 388
7800 AJ Emmen
The Netherlands

P. +31 591 631 700
F. +31 591 630 035
E. info@sparkholland.com
W. www.sparkholland.com

Head Office:
P. de Keyserstraat 8
7825 VE Emmen
The Netherlands





LC TROUBLESHOOTING

LC Method Scaling, Part II: Gradient Separations

If scaling isocratic separations is so simple, why is gradient scaling so confusing?

In last month's "LC Troubleshooting" installment (1), we looked at how to scale isocratic separations when the column size or packing particle size is changed. The process is quite simple. First, find a column with approximately the same plate number, then adjust the flow rate to give an acceptable pressure. The most common problems that result from mistakes in this scaling process give somewhat lower resolution than is expected or higher pressures. With gradients, unexpected consequences may occur from the changes that may be relatively unimportant in isocratic methods. In this month's discussion, we turn our attention to the proper scaling of gradient methods.

Resolution and Plate Number

Last month (1) we looked at the fundamental resolution equation for isocratic conditions:

$$R_s = \frac{1}{4} N^{0.5} (\alpha - 1) (k / [1 + k]) \quad [1]$$

where R_s is resolution, N is the column plate number, α is the separation factor, and k is the retention factor. A similar equation can be stated for gradient separations:

$$R_s = \frac{1}{4} N^{*0.5} (\alpha^* - 1) (k^* / [1 + k^*]) \quad [2]$$

where N^* is the effective plate number under gradient conditions, α^* is the gradient separation factor, and k^* is the gradient retention factor. As with isocratic conditions, we must be careful to keep from changing the chemistry of the system by keeping the same brand and series of column packing, the same mobile phase, and the same

column temperature. We'll see below that we have some additional factors to be careful of with gradients. If we keep these things constant, k^* and α^* (the ratio of k^* values for two adjacent peaks) should remain constant. If α^* is unchanged, we will obtain the same resolution if we keep the same column plate number.

The column plate number cannot be measured easily under gradient conditions, so we measure it under isocratic conditions. Because the plate number is a characteristic of the column, the use of isocratic conditions is not a problem. We use the same approach we used with isocratic conditions to select a column with an equivalent plate number so that we maintain the same resolution with the scaled method. We saw that we could obtain approximately the same plate number if we kept the column length-to-particle-diameter ratio constant within a range of +50% to -25%. Thus, we can determine the desired column length from

$$L_2 = L_1 d_{p2} / d_{p1} \quad [3]$$

where L_1 and L_2 are the column lengths, and d_{p1} and d_{p2} are the particle diameters of the original and new column, respectively. So far, nothing is different between the isocratic and gradient scaling process.

Gradient k^* Is the Key

With isocratic separation, we saw that the flow rate had no influence on the retention factor, k , and only a minor influence on the plate number, N . So, although we often start the isocratic scaling process by keeping the linear velocity of the mobile phase the same,


The Only Universal LC Detector Your Lab Will Ever Need

See What Other Detectors Are Missing

Charged aerosol detection is a revolutionary technology that will change the way you view every sample. This technique delivers consistent analyte response independent of chemical characteristics over a wide dynamic range while providing sensitivity at sub-nanogram levels.

Improve Inter-Analyte Response

An analyte's response to charged aerosol detection does not depend on optical properties, light scattering or the ability to ionize. Chromophores, radiolabels, ionizable moieties, or chemical derivatization are not essential for detection. Charged aerosol detection is a mass-sensitive technique that measures any non-volatile and many semi-volatile analytes. Variance in inter-analyte relative response is minimal, whether analyzing small molecules or proteins. And, this technique is gradient compatible.



Six pharmaceutical agents with an excipient (citric acid) were fully resolved using gradient reversed-phase HPLC and their responses measured first by UV detection and then by charged aerosol detection. As can be seen, UV detection significantly underestimates the levels of most analytes.

Unbiased Universal Detection

Charged aerosol detection has the flexibility to be used for a broad range of analytes in many different matrices, opening new opportunities for broad discovery and enhanced routine analysis.

- Drugs, impurities, and contaminants
- Biomolecules
- Foods and beverages
- Natural products, supplements, and botanicals
- Specialty chemicals
- Surfactants and polymers

Easy Integration With Any LC System

The new Thermo Scientific™ Dionex™ Corona™ Veo™ detector is designed to integrate into any HPLC/UHPLC system. When combined with a UV, diode array, or mass spectrometer, it provides an orthogonal and complementary detection solution, making it the ideal detector for any laboratory.



Reliable Results Without Intricate Optimization

The Corona Veo charged aerosol detector delivers sensitive, universal response through a simple yet flexible design, perfectly matched for applications with capillary, microbore, and analytical scale columns.

Download an application guide
or watch a video and see how
charged aerosol detection works:
thermoscientific.com/Veo



from a practical standpoint the flow rate is adjusted to maintain an acceptable pressure. With gradients, however, a change in the flow rate without some other compensating change can generate problems. This is because k^* is affected by flow rate:

$$k^* = (t_G F) / (1.15 V_m \Delta \Phi S) \quad [4]$$

where t_G is the gradient time (minutes), F is the flow rate (milliliters per minute), V_m is the column volume (milliliters), $\Delta \Phi$ is the gradient range (for example, 5–95% = 0.9), and S is a constant related to the molecular weight and is characteristic of each analyte. For the present discussion on scaling separations, we will not be changing the gradient range or the analytes, so $\Delta \Phi$ and S will be constant and can be dropped for a simpler relationship:

$$k^* \approx (t_G F) / V_m \quad [5]$$

With isocratic separations, the separation factor, α , is an expression of the peak spacing of two peaks with k values of k_1 and k_2 :

$$\alpha = k_2 / k_1 \quad [6]$$

As long as we do not change the chemistry of the system (mobile phase composition or column chemistry) or the column temperature, α will remain constant, as will the peak spacing. Isocratic k is calculated as follows:

$$k = (t_R - t_0) / t_0 \quad [7]$$

where t_R is the retention time and t_0 is the column dead time. If t_R is changed by nonchemical or temperature means, t_0 will change in parallel, so k stays constant. This is why we can change the column length, the column diameter, or the flow rate and not affect k in isocratic separation. If k is constant, α will also remain constant. On the other hand, if we do something that changes k , we expect that α will change because it is rare for a variable that changes k for one compound to change k in exactly the same manner for all the other compounds in the separation. Thus, a change in the mobile phase or column temperature usually results in

a change in peak spacing — in fact, we count on this effect as a tool for moving peaks relative to each other during method development.

A similar requirement of keeping α^* constant holds for gradient separations, where α^* is expressed as follows:

$$\alpha^* = k_2^* / k_1^* \quad [8]$$

As with isocratic separation, chemical or temperature changes are not likely to change k_1^* and k_2^* in the same way, so these variables must be kept constant with gradients too. But with gradients, the added complication is that the variables on the right-hand side of equation 5 must also be kept constant, or they should be changed in a manner that allows α^* to stay constant. Next, let's consider how this requirement influences the scaling of gradient methods.

Scaling a Simple Gradient

We'll start with a simple linear gradient of 5–95% acetonitrile over 20 min at 1 mL/min with a 150 mm \times 4.6 mm, 5- μ m particle diameter (d_p) column. Let's say that we want to get the same separation under ultrahigh-pressure liquid chromatography (UHPLC) conditions on a 1.8- μ m d_p column with an internal diameter of 2.1 mm.

The first step is to find a column with approximately the same plate number. We can use equation 3 to help us. The new column length (L_2) should be (150 mm \times 1.8 μ m) / 5 μ m = 54 mm long. Because 54-mm columns aren't available, we can use any length between 41 mm (-25%) and 81 mm (+50%) and be within our -25% to +50% guidelines. I would probably choose a 50 mm \times 2.1 mm column in the present case.

As with isocratic separations, with gradients it usually is a good idea to start with a flow rate that gives us the same linear velocity. We use the same technique we discussed last month:

$$F_2 = F_1 (d_{c2} / d_{c1})^2 \quad [9]$$

where F is the flow rate and d_c is the column internal diameter for columns 1 and 2. In the present case, the new flow rate should be 1 mL/min \times (2.1

mm / 4.6 mm)² = 0.208 mL/min; we'll round this to 0.2 mL/min for convenience. So, our new method will run on a 50 mm \times 2.1 mm, 1.8- μ m d_p column at 0.2 mL/min.

Let's check to see if we are done. Because we want k^* to remain constant, we can restate equation 5 as follows:

$$(t_{G1} F_1) / V_{m1} = (t_{G2} F_2) / V_{m2} \quad [10]$$

We can calculate each half of equation 10 for the proposed conditions and see how close we are. For the column volume, we'll use the approximation of ($L d_c^2$), because these are the only two variables that change between the columns. For the original conditions, (20 \times 1) / (150 \times 4.6²) = 0.00630, and for the proposed conditions, (20 \times 0.2) / (50 \times 2.1²) = 0.0181. The value for the new conditions is approximately three times too large. Because we've already defined the column size and the flow rate, the only variable we can adjust to fix this problem is the gradient time. If we reduce the gradient time to 7 min, the calculated value drops to 0.0635, which is close enough. So the gradient equivalent to the original conditions is a 7-min gradient at 0.2 mL/min on a 50 mm \times 2.1 mm, 1.8- μ m d_p column. This means that the switch from the high performance liquid chromatography (HPLC) to the UHPLC column reduces the run time by approximately threefold (20 min to 7 min).

But if we're running under UHPLC conditions, we can operate at higher pressures and because we're using sub-2- μ m particles, flow rate will not affect N , so a higher flow rate may help to further speed the separation. We use the same equation to calculate the system pressure as we did last month for isocratic conditions:

$$P_2 = P_1 (L_2 / L_1) (d_{c1} / d_{c2})^2 (d_{p1} / d_{p2})^2 (F_2 / F_1) \quad [11]$$

Where P_1 and P_2 are the original and new pressures, respectively; all the other variables were defined above. Let's assume that the original method generated a 2000 psi (140 bar) back pressure. For the proposed conditions, $P_2 = 2000 \times (50/150) \times (4.6/2.1)^2$



No Limits...No Compromises

A New Standard of UHPLC/HPLC Excellence

Utilizing patented technologies, and driven by customer input, Shimadzu has developed the **Nexera X2 Series** to deliver superior, real-world performance across a wider application range, while delivering unparalleled flexibility and reliability in UHPLC/HPLC analyses. With newly added features, such as **i-PDeA** and **i-DReC**, Nexera X2 will improve your productivity while assuring data quality. Nexera X2 is the next milestone in the evolution of liquid chromatography.

Learn more about Shimadzu's Nexera X2.
Call (800) 477-1227 or visit us online at
www.ssi.shimadzu.com/X2

Order consumables and accessories on-line at <http://store.shimadzu.com>
Shimadzu Scientific Instruments Inc.
7102 Riverwood Dr., Columbia, MD 21046, USA

Shimadzu's Nexera X2 HPLC/UHPLC features:

- **New PDA detector:** offers superior sensitivity and resolution, achieving a 0.4×10^{-5} AU noise level for genuine UHPLC analysis
- **Automated solvent blending and method scouting:** ideal for method development
- **Intelligent Peak Deconvolution Analysis (i-PDeA) function:** unique software solution to separate and accurately integrate co-eluting peaks
- **Intelligent Dynamic Range Extension Calculator (i-DReC):** extends the dynamic range so low concentration (trace or impurity compounds) and very high concentration samples can be accurately analyzed in a single injection
- **Ultrahigh-throughput and clean autosampler:** perfect front-end systems for any LCMS applications

Nexera X2
Ultra High Performance Liquid Chromatograph

$\times (5/1.8)^2 \times (0.2/1) \approx 4935$ psi. For UHPLC, we probably can tolerate twice this pressure, and the easiest way to do this is to double the flow rate so that a new pressure of ~ 9875 psi would be expected. Wait! We've changed the flow rate, so equation 5 tells us that this will change k^* , and we may find a change in peak spacing — not something we want. So we need to make another compensating change using equation 5 as a guideline. If we double F , we can reduce t_G by the same factor of two and k^* will stay constant. Now our method would be a 3.5-min gradient at 0.4 mL/min on a 50 mm \times 2.1 mm, 1.8- μ m d_p column. This would reduce the original run time by approximately sixfold (20 min to 3.5 min).

What About Segmented Gradients?

The example above used a simple, one-segment linear gradient. Let's add a bit of complication and see how

a two-segment gradient is handled. We'll use a 5–50% gradient in 10 min followed by a 50–95% gradient in 5 min, the other conditions (column and flow rate) are the same as before. Let's scale this to the same 2.1-mm i.d., 1.8- μ m d_p column. The same column, flow rate, and pressure scaling results will be obtained: scale to a 50 mm \times 2.1 mm column operated at 0.2 mL/min resulting in a back pressure of ~ 4935 psi.

When it comes to scaling the remaining gradient conditions, we proceed in the same manner as with a single-segment gradient, but we need to treat each segment as a separate gradient. Thus, we have two gradients to consider: 5–50%/10 min and 50–95%/5 min. Because we're not changing the gradient range, we can scale the gradient times so that equation 10 is balanced. For the first segment, the original conditions give $(10 \times 1)/(150 \times 4.6^2) = 0.00315$ and the proposed conditions give $(10 \times 0.2)/$

$(50 \times 2.1^2) = 0.00907$, approximately three times too large. A change in the gradient time to 3.5 min gives a value of 0.00317. For the second segment, our original gives $(5 \times 1)/(150 \times 4.6^2) = 0.00158$. We would expect the same gradient time ratio adjustment to hold for the second segment as for the first, so we can use it in our trial calculation. So we'll use $(3.5/10) \times 5$ min = 1.75 min, rounded to 1.7 min for convenience: $(1.7 \times 0.2)/(50 \times 2.1^2) = 0.00154$, which is close enough. Now our new method will be 5–50% in 3.5 min followed by 50–95% in 1.7 min at 0.2 mL/min on a 50 mm \times 2.1 mm, 1.8- μ m d_p column.

If we want to increase the flow rate to accommodate the higher pressure capability of a UHPLC system, we use the same technique as discussed above with equation 11, but remember to adjust the gradient time of each segment by a factor of two using equation 5 if you change the flow rate by a factor of two.

HPLC 2014



41st International Symposium on
High Performance Liquid Phase
Separations and Related Techniques

SYMPOSIUM CHAIR: Professor J. Michael Ramsey

SYMPOSIUM ORGANIZER: Professor Edward Yeung

SYMPOSIUM/EXHIBIT MANAGER: Ms. Janet Cunningham, Barr Enterprises

CONTACT

Ms. Janet Cunningham
Symposium/Exhibit Manager
Barr Enterprises

www.Linkedin.com/in/BarrEnterprises

@HPLC2014

janetbarr@aol.com

www.HPLC2014.org

Don't miss the largest global, multidisciplinary gathering of separation scientists at HPLC 2014, the conference of the year that covers all areas of liquid-phase separation science and related techniques.

Join us for an exciting program where experts from around the globe will gather to share their most significant, successful accomplishments and challenges facing their organizations.

If you want to learn about new methods, preview emerging technologies, ask questions to our technology experts, network to solve practical and economical issues, or hear about the latest innovations in all types of liquid-phase separations and related techniques, then we encourage you to plan to attend HPLC 2014 New Orleans.



May 11-15, 2014 • Hilton New Orleans Riverside • New Orleans, Louisiana, USA

Don't Forget to Adjust the Injection Volume

One thing I didn't mention in last month's discussion of isocratic scaling is the potential need to adjust the injection volume when the column size is changed. The same rules hold for both isocratic and gradient conditions. As a general rule, the injection volume should be adjusted in proportion to the column volume. We can use the following relationship as a guideline:

$$V_{inj2} = V_{inj1}(L_2 d_{c2}^2)/(L_1 d_{c1}^2) \quad [12]$$

Where V_{inj1} and V_{inj2} are the injection volumes of the original and new conditions, respectively. If our original method used a 25- μ L injection, the new method should use an injection of $(25 \times 50 \times 2.1^2)/(150 \times 4.6^2) = 1.75 \mu\text{L}$, which we can round to 2 μL .

Often the injection volume for the original method may not have been optimized, so although equation 12 gives the proper scaling, I suggest using an empirical test to determine the injection volume of the new method. I would try the calculated volume plus additional injections at twice and half the recommended volume. In the present case, I would try injections of 1, 2, and 4 (or 5) μL and observe the chromatogram. If there is no unacceptable degradation of the chromatogram, especially in the first part of the run, it is likely that the injection volume can be tolerated. I like to have a safety margin in my conditions, so the injection volume should not be so large that it is right on the edge of failure. This can be achieved by allowing a safety margin of 50–100%. For example, if in the above example, 5 μL still looked OK, try 8 or 10 μL . If these injection volumes are OK, 5 μL is safe.

Further Cautions

As you can see from the discussion above, when the number of segments in the gradient increases, the work in calculating new conditions increases because each gradient segment must be scaled individually. What happens when there is a curved (convex or concave) gradient? I recommend strongly against using this type of gradient, because they are hard to reproduce from one instrument to another, and very difficult to scale. I don't know of any way to scale such gradients, because they are essentially composed of an infinite number of very small segments with different slopes. So, if you have an existing curved gradient that needs to be scaled, use the single-segment rules and hope that the new gradient works.

Those of you who are familiar with gradients will notice that I ignored the instrument dwell volume (gradient delay volume) in this discussion. The new dwell volume requirements can be calculated if the dwell time (the time to wash out the dwell volume) is considered as a gradient segment. The complication is that the dwell volume is seldom adjustable. Usually when you are moving a method from a conventional column to a small-volume or UHPLC column, you also will be moving from a conventional LC system to a UHPLC or a newer LC system, either of which is likely to have a reduced dwell volume. Dwell volume differences are most likely to affect the peak spacing for peaks eluted early in the gradient, so watch for such changes when the method is scaled. You may need to adjust the initial

conditions to compensate for such changes. Dwell volume, of course, is not important in isocratic methods.

Finally, as with the isocratic examples last month, gradient calculations are tedious. You can simplify the process by using one of the on-line calculators to do the work for you. Some of these also calculate scaling of the injection volume and take dwell volume into account. Search the internet for "HPLC method transfer calculator" and you'll find several to choose from. I suggest that you try several and settle on the one that seems like the easiest for you.

References

- (1) J.W. Dolan, *LCGC North Am.* 32(2), 98–102 (2014).

John W. Dolan


"LC Troubleshooting" Editor John Dolan has been writing "LC Troubleshooting" for LCGC for more than 30 years. One of the industry's most respected professionals, John is currently the Vice President of and a principal instructor for LC Resources in Walnut Creek, California. He is also a member of LCGC's editorial advisory board. Direct correspondence about this column via e-mail to John.Dolan@LCResources.com



For more information on this topic,
please visit www.chromatographyonline.com/dolan

ACE® UltraCore™

All the proven advantages of ACE UHPLC / HPLC columns...
...now available with solid-core particles



- Ultra-inert solid-core UHPLC / HPLC columns
- High efficiency 2.5 μm and 5 μm solid-core particles
- SuperC18 and SuperPhenylHexyl phases with extended pH stability



Contact your local distributor
or visit www.ace-hplc.com

ACE, UltraCore, SuperC18 and SuperPhenylHexyl are trademarks of Advanced Chromatography Technologies Ltd



GC CONNECTIONS

The Storage and Use of Gases for Gas Chromatography

This month's column discusses the safe storage and use of calibration, carrier, and detector gas cylinders for small or large organizations.

Carrier and detector gases, when not generated directly on-site, are transported to laboratory facilities in gas cylinders and stored there until they are used. A recent "GC Connections" installment (1) explored the sources, both natural and synthetic, of carrier and detector gases for gas chromatography (GC). When put into service, secured gas cylinders or laboratory generators are connected to pressure- and flow-regulating devices that establish appropriate flow and pressure conditions for delivery to one or more gas chromatographs. The delivery system may be a simple length of tubing or it may be a more complex manifolded arrangement with intermediate regulation and valving. Ideally, the gas stream passes through a final purification stage positioned close to each chromatograph.

Much of the equipment associated with gas supplies in the laboratory ensures the purity and control of carrier and detector gases; the rest provides necessary safety measures for day-to-day encounters with flammable or asphyxiant gases that come in heavy high-pressure containers. Both the correct equipment and the procedures associated with gas deployment are critical for achieving the safest work environment and the best possible results.

In the United States, both the federal and state government authorities promulgate regulations and guidelines for the safe handling and use of gases according to their delivery system and chemical nature. In 29 *CFR* 1910.101, the United States Occupational Health and Safety Administration (OSHA) incorporates standards issued by the Compressed Gas Association (CGA) (2). Several organizations have made their policies and procedures for compressed gas handling available on-line; among these, the State University of New

York (SUNY) campus at Stony Brook Compressed Gas Safety Guide (3) is one that is easily accessible and useful for typical chromatography laboratory scenarios.

Industrial gas consumers — including both companies and their employees — are responsible for the safe storage and use of gas cylinders and their contents from the moment the tanks come off a delivery truck until the gas supplier collects the empties for return. During this period of responsibility, certain procedures and equipment are used to mitigate the hazards and get the best use out of the cylinders' contents. In addition, companies should enact specific training programs to qualify employees to handle compressed gases safely and effectively. The contents of this article are intended as general professional advice and are not meant to reduce or replace any requirements, procedures, or regulations imposed by companies or their local authorities.

Gas Cylinders

Gas cylinders are ubiquitous in chromatography laboratories. Their presence in nearly all chromatography work environments, combined with the rarity of accidental injury or property damage, engenders a sometimes too-relaxed attitude toward them on the part of laboratory workers and managers alike. Other hazards, such as burns from hot GC inlets or punctures suffered while cutting fused-silica columns, are much more common and can be mitigated by following easily understood and simple procedures. It is safe to say that the gas cylinders in a chromatography laboratory present an array of hazards that, short of the presence of toxic chemicals, comprise the greatest personal and property safety challenges chromatographers are likely to encounter.

evolve into GC 2.0

Unique and innovative breakthroughs featured on the **Thermo Scientific™ TRACE™ 1300 Series Gas Chromatograph** elevate performance and empower the user. Developed around customer-driven innovations, including user-installable instant connect injectors and detectors, unique and patented solutions to conserve helium carrier gas while running analyses, and miniaturized electronic components, **the TRACE 1300 Series GC is an extremely fast, easy to use, compact GC**, delivering incredibly high lab productivity at a much reduced cost of ownership.

powerful breakthroughs

- see the best choice in gas chromatography at thermoscientific.com/trace1300



TRACE 1300 Series GC
Transcend traditional GC

- Instant connect, user-installable injectors and detectors
- Ultimate productivity in routine analysis
- Fast, high-sensitivity detectors



Instant connect module installation by user



Thermo Scientific™ Dionex™ Chromeleon™ Chromatography Data System software for complete control and compatibility



Thermo Scientific Instant Connect Helium Saver technology eliminates helium supply shortages and saves money while keeping GC and GC-MS methods intact

Cylinder Hazards

Gas cylinder users do not always follow safe procedures. I often see cylinders that are not properly restrained or lack valve-protection caps that should be installed when the cylinders are not in use. Sometimes I observe incorrect transfer of cylinders from one location to another. For example, years ago at a university a gas delivery worker tilted two cylinders 15° from vertical, one with each gloved hand, and proceeded to roll them on their bottom edges down a hallway and into a laboratory. The cylinders did not fall, but I wonder how many times they had fallen during his tenure in that job and what kind of damage had been caused.

And how about moving tanks with a regulator attached and the main valve open, or without a protective valve cap? Not a good idea, but I've seen it done more than once. So, why not handle and use compressed gas cylinders without following the prescribed safe procedures? There are lots of reasons not to do so, and they can be grouped into three categories: physical hazards, energy hazards, and chemical hazards.

Physical Hazards

Gas cylinders are massive. A 49-L capacity steel cylinder, such as what is commonly used for helium or nitrogen carrier gas, weighs about 63 kg (138 lb) and is physically unstable in its normal vertical position. Such a mass of steel can do considerable damage to limbs and toes that happen to be in the way as it falls.

Energy Hazards

In addition to being bulky, massive, and unstable, compressed gas cylinders constrain a considerable amount of potential energy as pressurized gas. Should that pressure be released suddenly, a cylinder can become a spinning projectile that can penetrate walls and potentially inflict deadly harm. In an earlier "GC Connections" installment (4), I calculated that a full helium cylinder could attain a velocity of more than 108 km/h (67 mph) in a few seconds if the gas were released along the cylinder's main axis. The result depends on how quickly the gas is released, but by an estimate that is enough energy to do a lot of damage.

Chemical Hazards

In the case of hydrogen, a gas cylinder represents considerable combustible chemical energy as well. Hydrogen concentrations can reach the lower explosive limit (LEL) of only 4% in air in any enclosed space into which sufficient quantities are released.

Cylinders that contain gases other than air also pose an asphyxiation hazard that results from the potential displacement of oxygen should the gas be released more rapidly than ambient air can be replaced by ventilation.

Some laboratories also use cryogenic liquefied gases such as liquid nitrogen or carbon dioxide. The proper use and storage of cryogenic liquids is a separate topic about which the interested reader can learn more from the references at the end of this article. Specific personnel protective equipment such as face masks, gloves, and boots should always be used when handling cryogenic liquids.

Built-In Protection

Gas cylinders for the laboratory include some built-in safety measures. They are protected against over-pressurization caused by heating or inadvertent over-filling by a burst or rupture disk in the cylinder valve stem. Safety burst disks will open at pressures greater than the normal operating pressure for the cylinder, but less than the cylinder's test pressure. For a laboratory carrier-gas tank, the disk will open when the tank pressure exceeds about 225 MPa (3200 psi). This pressure is greater than what would be attained should a completely full cylinder be stored at up to 60 °C; the maximum recommended storage temperature for nonreactive gases is around 50 °C. A ruptured disk will release the contents of the cylinder safely (and loudly). A so-affected cylinder must be returned to the supplier for inspection and refurbishing.

Managing Cylinders

Now let's examine how to manage and handle gas cylinders safely. A cylinder passes through three stages at a laboratory site: receipt and storage, in-use, and ready to return. These stages usually represent three separate locations within a facility; the first and last have nearly identical requirements, while the in-use stage is different.

Regardless of the stage of a cylinder's lifetime, store inert laboratory gases separately from reactive gases, including hydrogen fuel gas or any other gases used for other purposes. Hydrogen has some special venting and storage requirements that are detailed in the standards; keeping different types of gases separate helps avoid mix-ups.

MARCH 24-25, 2014
SHERATON HOTEL
 SILVER SPRING, MARYLAND USA



BIOASSAYS
 SCIENTIFIC APPROACHES AND REGULATORY STRATEGIES 2014



TOPICS THAT WILL BE DISCUSSED:

- Bioassays Lessons Learned Case Studies
- Bioassays to Support Biopharmaceutical Development
- Bioassay Controls & Control Strategies
- Vendor Showcase



SCAN
to browse
CASSS
for program
updates at casss.org



CASSS
INTERNATIONAL SEPARATION SCIENCE SOCIETY
www.casss.org

Receipt and Storage

Receiving personnel should check the labeled contents of all cylinders upon delivery, as well as the overall appearance and condition of the cylinders. Any exceptions should be brought to the attention of the supplier as soon as possible. Suppliers don't deliver damaged or incorrect cylinders intentionally, but a quick check is in order. Ideally, the date of the most recent hydrostatic pressure test should be verified; however, this information can be difficult to ascertain from the markings on a cylinder. Cylinder pressure can be verified if the valve assembly includes an indicating regulator; otherwise it's better to wait until a cylinder is put into service before checking the fill pressure.

Gas cylinder contents generally are color-coded on the cylinder, but the color scheme is not consistent from supplier to supplier. Thus, it is not possible to rely on the color of a cylinder for identification. An attached label or tag will give the necessary information. If there is no tag or label, then the cylinder must be returned immediately to the supplier — or another supplier if the original supplier is unknown — for proper disposition.

Gas cylinders must be moved and stored properly from receipt until they are returned to the supplier. Outdoor and indoor storage are both possible. Indoor storage is more convenient if space is available, while outdoor storage should provide protection from weather and dirt. Cylinders must be kept dry in all cases to avoid corrosion.

Suitable cylinder storage locations include any area with limited access where they can be properly restrained. For small organizations the laboratory itself is often the best spot. Storage in a public hallway is not desirable, of course, while too-limited access might cause a problem should quick access be required in an emergency. Each location where gas cylinders are used is unique and, if at all possible, the question of proper cylinder locations should be considered carefully.

Often, and especially in situations where a larger number of cylinders are used, there may be more than one location for cylinder storage. For example, some companies store full cylinders in a location near a loading dock while keeping in-use cylinders in the laboratory or in a shared area that feeds multiple locations.

In-Use

Cylinders are classified as in-use if a regulator is attached or the protective valve cap is removed, whether or not gas is flowing at any particular time. With the protective cap removed or with a regulator attached, the high-pressure cylinder valve is exposed and represents an increased hazard. A cylinder that falls on its side with the cap or protective collar in place is unlikely to decompress explosively, whereas one with an exposed valve is much more likely to do so. Therefore, the in-use category requires very careful attention to cylinder restraints.


Ready to Return

Used cylinders ready for return should be kept in a dedicated empty-cylinder area that is convenient for their removal from the building, or at least the empties should be marked clearly as such. Once empty or otherwise no longer needed, cylinders fall back into a category with the same requirements as full cylinders waiting to be put in use — the protective valve cap must be installed; otherwise the valve could be damaged. Even nonobvious damage to the valve can



COSMOCORE

Increased loading capacity
Excellent pH stability (1-10)



nacalai tesque
The quality for certainty.



www.nacalaiusa.com

be a problem for the supplier when the cylinder is refilled.

If possible, gas cylinders should not be emptied completely before they are returned, since doing so may expose the inside of the cylinder to ambient air and moisture. The worst case is for an empty cylinder to be stored outside with an open valve. It will be difficult for the supplier to completely purge and clean the cylinder before the next filling. A residual pressure of 175–700 kPa (about 25–100 psi) is suitable.

Moving Gas Cylinders

The image of two large cylinders being rolled by hand on their edges down a hallway remains with me more than 30 years after I witnessed it. The right way to move cylinders is of course with a cylinder cart or hand-truck. These are available in single- and dual-cylinder models with a variety of configurations. The better designs include static-dissipating wheels, which are much appreciated when moving hydrogen cylinders. Never move a cylinder with the regulator attached or the valve

cover missing. And never use a gas cart as a permanent cylinder stand. Some tilting and rolling is required to maneuver a cylinder onto the cart, so caution is appropriate. Wear gloves and protective goggles and footwear.

If a cylinder starts to fall from upright toward prone, do not attempt to catch it. Jump back quickly out of the way. Return a fallen or otherwise physically stressed cylinder to the supplier immediately. A fall can introduce weaknesses in the cylinder that may make it potentially unsafe for continued use.

How to Restrain Gas Cylinders


Proper restraint of gas cylinders is extremely important for their safe use. When set into position, cylinders must be prevented from tipping over, especially when in use with regulators attached. Gas cylinder restraints provide protection from accidental tippage that could compromise attached gas lines or, in the extreme case, cause a regulator or the tank valve to break or detach.

Several types of restraining systems are available, including single- or dual-tank straps that attach to a fixed object or a wall; individual free-standing cylinder floor stands; small multicylinder corrals with a restraining chain; and completely enclosed and lockable cylinder cages.

Do not place cylinders on their sides: moving them down to the floor or back to the upright position is a lifting hazard at the least, and the movement greatly increases the chances for dropping the cylinder.

Tank Straps

These are available in several forms. They all feature an attachable strap or brace that encircles one or two cylinders approximately 1 m above the floor. The chain, woven belt, or metal strap is attached with a snap, buckle, or bolt so that the cylinder is prevented from falling over or moving out of position. Tank straps should be affixed permanently to an immovable object such as a wall or a laboratory bench that is bolted to the wall or floor. Some straps feature a clamp that attaches temporarily to a benchtop and allows the strap to be relocated easily. Although this type does not provide strict immobility it may be sufficient in some locations, but only if the bench is secured.



**SEPARATE
PEPTIDES & PROTEINS**
faster
with higher peak capacity
at lower back pressure
using higher temperature

HALO[®] BIOCLASS
PEPTIDE ES-C18 • PEPTIDE ES-CN
PROTEIN C4 • PROTEIN ES-C18

The same Fused-Core[®] technology that led the revolution in core-shell columns is now available in HPLC and UHPLC columns specifically designed for the separation of peptides and proteins.

For more information, visit www.advanced-materials-tech.com

© HALO and Fused-Core are registered trademarks of Advanced Materials Technology, Inc.

Floor Stands

Floor stands are single-cylinder stands that clamp or strap onto a cylinder and provide a larger footprint than the cylinder itself. They provide improved stability but are not immune to tipping or movement, so they are recommended only for temporary cylinder placement.

Cylinder Corrals

A cylinder corral is affixed permanently to the floor or wall. The corral admits two or three cylinders across and two or three deep, and the cylinders are restrained by a chain that attaches onto several hooks welded to the corral. This arrangement is suitable for stored or in-use cylinders where single tank straps are not enough — for example, in a gas analysis laboratory where multiple calibration gas cylinders plus carrier gas, hydrogen, and air are required. A cylinder corral is the most flexible arrangement and requires the least amount of linear wall space for multiple cylinders.

Cage Enclosures

A cylinder cage or locker is a large box with heavy mesh walls and usually a floor as well. The door is firmly secured with a lockable bolt. Cylinder cabinets are similar, but they have solid walls. A cylinder cage holds up to 16 or more laboratory gas tanks, but only for storage purposes as it is not a good practice to snake gas lines through the cage mesh. It is the most secure method of cylinder storage. Cages or cabinets are also the most expensive storage method, but they are well worth the investment in safety and security for organizations that have to manage a large number of cylinders.

Using Laboratory Gases

After verifying gas cylinders and their content upon receipt, safely moving and securing them in place, and organizing them by their current use, the next step is to attach the appropriate gas regulators, tubing, filters, and fittings to bring the high-purity gas supplies up to one or more gas chromatographs. The next "GC Connections" installment will discuss these topics and more, as we continue to follow the multiple gas paths through a GC system.

References

- (1) J.V. Hinshaw, *LCGC North Am.* **32**(1), 30–35 (2014).
- (2) Compressed Gases, in 29 *CFR* 1910.101 (U.S. Government Printing Office, Washington, D.C.). Available at https://www.osha.gov/pls/oshaweb/owadisp.show_document?p_table=standards&p_id=9747.
- (3) "Compressed Gas Safety Guide," at <http://www.stonybrook.edu/facilities/ehs/occupational/cg.shtml> (SUNY Stony Brook University, January 2014).
- (4) J.V. Hinshaw, *LCGC North Am.* **20**(6), 532–536 (2002).

John V. Hinshaw

"GC Connections"
editor John V. Hinshaw
is a Senior Scientist at
Serveron Corporation in
Beaverton, Oregon, and
a member of LCGC's editorial advisory board.
Direct correspondence about this column to
the author via e-mail: lccedit@lccgmag.com.



For more information on this topic,
please visit
www.chromatographyonline.com/hinshaw

O.I. Analytical
a xylem brand

INNOVATIVE
SAMPLE
PROCESSING FOR
VOC ANALYSIS

The 4100 Water / Soil Sample Processor efficiently processes up to 100 samples for purge-and-trap analysis of VOCs. USEPA -approved methods for drinking water, wastewater, or soil samples can be run with the 4100 connected to one or two purge-and-trap instruments.

The innovative VOA Constrictor™ cylindrical gripper lifts and transports 40-mL VOA vials with unsurpassed reliability.

www.oico.com / 800.653.1711

xylem
Let's Solve Water

Analysis of Amphetamines and Synthetic Cathinones in Hair Samples Using LC–Tandem Mass Spectrometry and Solid-Phase Extraction



In this study, the procedure for analyzing amphetamines and synthetic cathinones (also known as bath salts drugs) in hair samples using a mixed-mode solid-phase extraction (SPE) is described. Samples of hair were digested with a dilute solution of base (containing internal standards), neutralized, and diluted with an aqueous phosphate buffer (pH 6). Each sample was applied to a conditioned SPE column, after which the sorbent was rinsed with deionized water, acetic acid, and methanol. After drying, the analytes were eluted and collected from the SPE column with 3 mL of an elution solvent consisting of methylene chloride–isopropanol–ammonium hydroxide. To the eluate, 200 μ L of mobile phase was added and the samples were evaporated to the mobile phase for analysis by liquid chromatography–tandem mass spectrometry (LC–MS–MS). Chromatography was performed in gradient mode using a C18 column and a mobile phase consisting of acetonitrile and 0.1% aqueous formic acid. The total run time for each analysis was under 5 min.

Amphetamine (α -methylphenethylamine) (Figure 1) belongs to a class of compounds known as *sympathomimetic amines* (1). This class of drugs includes the illicit drugs methamphetamine, methylenedioxyamphetamine (MDA), and methylenedioxy-methamphetamine (MDMA), as well as ephedrine and pseudoephedrine, which can be found in over-the-counter medications. Amphetamine is administered as a prescription medication (for example, Adderall [Shire Pharmaceuticals]) for treating medical issues such as narcolepsy, obesity, or hypotension (2), whereas methamphetamine, MDA, and MDMA are considered controlled substances — that is, pharmaceuticals with little or no medical use.

Synthetic cathinones are derived structurally from the parent compound (Figure 2) and have become noticeable in the scientific literature in recent times because of the fatalities arising from administration (3–6).

The drugs are commonly referred to as “bath salts” because they were originally packaged with names such as “Ivory Wave” and marketed as “not for human consumption” or “research chemicals.” These drugs are now scheduled in the same way as other controlled substances.

The popularity of amphetamines is because of their euphoria effect and ease of synthesis. Their use or abuse is generally verified by the analysis of biological samples, such as urine, blood, oral fluid, or hair. Of these samples, hair is a biological matrix that has been used as an alternative to urine or blood for drug testing because it allows noninvasive sampling and can document the use of the drugs over a longer period of time than blood or urine (7). In this study, amphetamine in the form of Adderall was determined in the hair of a subject along with several other amphetamines and a range of synthetic cathinones.

In terms of metabolism in the human system, amphetamine is inactivated during

Albert A. Elian*, **Jeffery Hackett[†]**, and **Kerrie T. Donovan***

*Massachusetts State Police Crime Laboratory, Sudbury, Massachusetts; [†]UCT Inc., Bristol, Pennsylvania. Direct correspondence to: jhackett@unitedchem.com




Figure 1: Structure of amphetamine.




Figure 2: Structure of cathinone.

this process, undergoing deamination to form phenylacetone, which is converted to benzoic acid and excreted in a conjugated form (2). A small amount of the parent is oxidized to norephedrine, which is also metabolized to the parahydroxylated forms

of this compound, all of which is pharmacologically active and are thought to contribute to the effects of the drug (8,9). The therapeutic, toxic, and fatal concentrations of amphetamines in samples such as blood and urine are well documented (2,10), but not so much for the synthetic cathinones because of the recent nature of their abuse. Most of the published studies on hair analysis have been performed in the area of workplace drug testing and drugs or driving cases (11,12), not postmortem studies.

Gas chromatography coupled to mass spectrometry (GC-MS) has been reported as a technique for quantifying amphetamines in hair (13). GC-MS analysis of the amphetamines requires derivatization using compounds such as heptafluorobutyric anhydride (HFAA) or pentafluoropropionic anhydride (PFAA) (14,15). Because amphetamine exists as a *d-l* isomeric pair, some laboratories have used chiral modification to separate the isomers in samples such as hair (16). The ratio of the isomeric forms may indicate whether or not the amphetamine has been taken legally. GC-MS analysis of the cathinones used similar fluoroacyl derivatives (3). Liquid chroma-

tography coupled to tandem mass spectrometry (LC-MS-MS) is gaining popularity for analysis in this matrix (17,18). The use of solid-phase extraction (SPE) described in this article uses the LC mobile phase as a keeper solvent for amphetamine, reducing its volatility. In previous methods the addition of methanolic hydrochloric acid or a solvent such as dimethylformamide has been reported, and the solvent was evaporated to dryness (19). The addition of the mobile phase presents the LC-MS-MS with a more amenable analytical solvent. SPE has been reported in the analysis of hair samples previously (20-22), but not using this type of keeper solution format.

Experimental

Chemicals and Reagents

Amphetamine, amphetamine-*d*₅, methamphetamine, methamphetamine-*d*₅, MDA, MDA-*d*₅, MDMA, MDMA-*d*₅, butylone, ethylone, flephedrone, mephedrone, mephedrone-*d*₅, methylone, methedrone, methcathinone (4-MEC), methylenedioxy-pyralerone (MDPV), and pyralerone were obtained from Lipomed as 1-mg/mL methanolic solutions. Acetonitrile, acetic

ANALYTICA
BOOTH A2.141

ASMS
BOOTH 54

TRAPS

Keep the Best • Flush the Rest

On-Line HPLC/UHPLC Trapping

For Small Molecule, Peptide and Protein Applications

PITCON
BOOTH 4048

MSACL
BOOTH 22




Table I: Tandem mass spectrometry conditions

Compound	Q1	Q3	t_R (min)	DP (V)	CE (V)	CXP (V)	CE (V)
Amphetamine (1)	136.1	91.0	1.30	51	10	23	16
Amphetamine (2)	136.1	65.0	1.30	51	10	13	8
Amphetamine- d_5 (2)	141.1	124.0	1.30	51	10	25	4
Amphetamine- d_5 (2)	141.1	93.0	1.30	51	10	21	6
Methamphetamine (1)	150.1	91.0	1.62	61	10	29	16
Methamphetamine (2)	150.1	119.1	1.62	61	10	17	8
Methamphetamine- d_5 (1)	155.2	92.1	1.62	46	10	57	10
Methamphetamine- d_5 (2)	155.2	121.2	1.62	46	10	25	8
MDA (1)	180.2	163.1	1.51	46	10	15	8
MDA (2)	180.2	105.1	1.51	46	10	33	8
MDA- d_5 (1)	185.2	168.1	1.51	46	10	33	6
MDA- d_5 (2)	185.2	110.1	1.51	66	10	55	14
MDMA (1)	194.2	163.1	1.73	76	10	19	12
MDMA (2)	194.2	105.1	1.73	76	12	33	18
MDMA- d_5 (1)	199.2	165.1	1.73	66	10	19	12
MDMA- d_5 (2)	199.2	106.8	1.73	66	10	19	12
Mephedrone- d_3 (1)	181.2	123.0	1.94	66	10	35	6
Mephedrone- d_3 (2)	181.2	163.0	1.94	46	10	17	10
Methedrone (1)	194.1	176.1	1.66	56	10	29	12
Methedrone (2)	194.1	161.1	1.66	56	10	39	10
Butylone (1)	222.1	174.2	1.56	35	5	22	26
Butylone (2)	222.1	204.2	1.56	35	5	22	34
Ethylone (1)	222.1	174.2	1.86	35	5	22	26
Ethylone (2)	222.1	174.2	1.86	35	5	22	34
Flephedrone (1)	182.1	164.2	1.22	37	5	18	19
Flephedrone (2)	182.1	149.1	1.22	37	5	18	30
Mephedrone (1)	178.1	145.1	1.94	36	4	18	28
Mephedrone (2)	178.1	160.1	1.94	36	4	18	18
Methylone (1)	208.1	160.1	1.15	66	10	27	10
Methylone (2)	208.1	132.1	1.15	66	10	19	12
4-MEC (1)	192.1	174.2	2.07	66	10	19	12
4-MEC (2)	192.1	144.1	2.07	66	10	41	10
MDPV (1)	276.2	135.1	2.32	101	10	37	12
MDPV (2)	276.2	126.1	2.32	101	10	31	16
Pyrvalerone (1)	246.2	105.1	2.55	56	10	33	188
Pyrvalerone (2)	246.2	175.2	2.55	56	10	59	16

DP = declustering potential; EP = exit potential; CXP = collision cell exit potential; CE = collision energy; t_R = retention time (in minutes).

acid (glacial), concentrated ammonium hydroxide solution (32% by volume), formic acid, isopropanol, methanol, and methylene chloride were obtained from Fisher Scientific. The SPE columns (CSDAU206)

were obtained from UCT Inc. Deionized (DI) water was laboratory grade, and it was generated in the Massachusetts State Police Crime Laboratory (MSPCL). All chemicals were of ACS grade.

Acetic acid was prepared as a 0.1 M solution by diluting glacial acetic acid (5.8 mL to 500 mL) and then increasing the volume to 1 L by adding DI water and mixing well. Formic acid was prepared as a 0.1%

BUILDING BETTER SCIENCE

ENJOY THE BENEFITS OF CHROMACADEMY

COMPLIMENTS OF AGILENT TECHNOLOGIES

Five years FREE access for all university students & staff at www.chromacademy.com/agilent

CHROMacademy is an intuitive, comprehensive e-learning and trouble-shooting platform with more than 3,000 pages of content for HPLC, GC, sample preparation, and hyphenated techniques. No other online resource offers separation scientists more live streaming events, a knowledge base, practical solutions, and new technologies in one easy to navigate website.

- Monthly webcasts increase your knowledge from the comfort and convenience of your desk
- Multi-media and video tutorials enhance your skills at your own pace
- Interactive experiments, lab simulations, and tools improve your productivity
- Applications, news, and feature articles from world leaders in analytical science
- Troubleshooting modules based on real-world laboratory experience, help you solve common problems
- Access to over 5,000 articles and application notes from LCGC
- ...and more.

In addition, CHROMacademy has an interactive discussion forum where visitors can get technical and application assistance from other users or provide advice to members who are looking for help.

All of this can be yours for **free***, thanks to Agilent Technologies. Simply log on to www.chromacademy.com/agilent, complete** and submit the form, and receive your **complimentary five year membership worth US \$1,475**.

* Five years free access to CHROMacademy only available to customers affiliated with an academic or research institution, conditions apply.


** A valid university e-mail address is required.

LC|GC's **CHROMacademy**



Stay connected with Agilent and CHROMacademy via Twitter, Facebook, LinkedIn and YouTube.

Join the conversation today!

 [@AgilentChem](https://twitter.com/AgilentChem) or [@AgilentLife](https://twitter.com/AgilentLife)

 [@CHROMacademy](https://twitter.com/CHROMacademy)

 linkedin.com/company/Agilent-Technologies

 YouTube.com/AgilentChem or YouTube.com/AgilentLife

 facebook.com/Agilent.Tech

 facebook.com/CHROMacademy

Agilent Products are for Research Use Only.
Not for use in diagnostic procedures.
Information, descriptions and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc. 2013
Published in USA, 5991-1198ENUC, 2013



Agilent Technologies

Table II: Precision results for control hair samples shown as percent recovery

Compound	2 ng/mg	8 ng/mg
Amphetamine	100 ± 5	102 ± 7
Methamphetamine	95 ± 5	101 ± 7
MDA	110 ± 15	103 ± 11
MDMA	105 ± 5	102 ± 10
Butylone	120 ± 20	106 ± 12
Ethylone	90 ± 15	108 ± 13
Flephedrone	85 ± 20	103 ± 8
Mephedrone	120 ± 15	106 ± 10
Methylone	115 ± 20	98 ± 10
Methedrone	106 ± 7	101 ± 13
4-MEC	90 ± 15	103 ± 12
MDPV	95 ± 10	97 ± 11
Pyralverone	115 ± 10	103 ± 10

(v/v) solution by adding 1 mL of the acid to 900 mL of DI water and diluting it to 1 L (mobile-phase A). Acetonitrile containing 0.1% formic acid (v/v) was prepared by adding 1 mL of formic acid to 900 mL of acetonitrile and diluting to 1 L (mobile-phase B). Aqueous sodium hydroxide was prepared as a 0.1 M solution by adding 4 g of the solid to 500 mL of DI water and dissolving before diluting to 1 L. Phosphate buffer (pH 6, 0.1 M) was purchased from Fisher Scientific as a ready-to-use solution.

Chromatographic Analysis

Analysis was performed using an API 3200 Q-Trap instrument supplied by Applied Biosystems. The chromatographic system consisted of a Shimadzu CBM 20 A controller, two Shimadzu LC 20 AD pumps (including a degasser), a Shimadzu SIL 20 AC autosampler, and a Shimadzu CTO AC oven (set at 10 °C). The instrument was fitted with a 50 mm × 2 mm, 5- μ m d_p Cadenza-C18 column from Imtakt USA (formerly Silvertone Sciences) that was attached to a 5 mm × 2 mm Cadenza-C18 guard column obtained from the same supplier. The LC column oven was maintained at 40 °C throughout the analyses. The injection volume was 10 μ L. The mobile phase was delivered at a flow rate of 0.5 mL/min. The mobile-phase gradient program started at 5% mobile-phase B, rose to 90% B in 4.0 min, and then returned to 5.0% B. The instrument was readied for reinjection after 5.1 min.

The mass spectrometry was performed using positive multiple reaction monitoring (MRM). The mass spectrometer conditions for each of the amphetamines and synthetic cathinones are shown in Table I. Tandem MS was performed using the following conditions: curtain gas setting, 15; collision gas setting, medium; ion spray voltage setting, 5000 V; temperature setting, 650 °C; ion source gas 1 setting, 50; ion source gas 2 setting, 50. The tandem mass spectrometer conditions are shown in Table I. The analytical data were collected using Analyst software version 1.5.2 supplied by Applied Biosystems.

Sample Preparation for Analysis

Calibrators and Controls

A solution of amphetamines and synthetic cathinones (amphetamine, methamphetamine, MDA, MDMA, butylone, ethylone, flephedrone, mephedrone, methylone, methedrone, methcathinone [4-MEC], methylenedioxypravalverone [MDPV], and pyralverone) was prepared at a concentration of 1 μ g/mL by the dilution of 10 μ L of stock solution with acetonitrile to 10 mL in a volumetric flask. A 1- μ g/mL solution of the internal standards (amphetamine- d_5 , methamphetamine- d_5 , MDMA- d_5 , MDA- d_5 , and mephedrone- d_3) was prepared by the diluting 100 μ L of the stock solution (100 μ g/mL) to 10 mL with acetonitrile in a volumetric flask. The choice of internal standard was based on the fact that deuterated analogues of amphetamines and synthetic cathinones would not only chromatograph in a very similar mode to the drugs themselves, but also would be extracted via SPE as efficiently as amphetamine and the synthetic cathinones and would not likely be observed in the samples under analysis.

Calibrators were prepared by the addition of 1.0, 2.5, 5.0, 10.0, 25.0, 50.0, and 100.0 ng of amphetamines and synthetic cathinones into 10-mg samples of drug-free hair samples. The hair samples had been previously decontaminated by washing two times with methylene chloride (10 mL) and two times with DI water (10 mL) before being air dried for 24 h. Each wash step was analyzed for the presence of drugs and found to be negative. Next, 50 ng of the internal standard was added to these samples. Control samples were prepared by the addition of 2 ng and 8 ng of amphetamine and synthetic cathinones to drug-free hair samples in addition to 50 ng

of the internal standard. All determinations were performed in duplicate. A negative control sample was prepared by the addition of only the internal standard (50 ng) to a sample of drug-free hair samples (10 mg). Calibrators, control samples, and test samples were treated in an identical mode with regard to sample extraction.

To assess the performance of the procedure, calibration curves were constructed twice daily over five consecutive days using the spiked controls, and from these data intraday and interday values were obtained.

Sample Extraction

A 1-mL volume of 0.1 M aqueous sodium hydroxide solution was added to each sample (calibrator, control, and test) in a clean glass tube (75 mm × 12 mm) that was then capped. The tubes and contents were digested for 30 min at room temperature. Then, 4 mL of 0.1 M phosphate buffer (pH 6) was added to the solution, which was mixed and centrifuged at 3000 rpm for 10 min before it was applied to the SPE column.

Solid-Phase Extraction

SPE columns were conditioned by the sequential addition of 3 mL of methanol, 3 mL of DI water, and 1 mL of 0.1 M phosphate buffer (pH 6). Each liquid was allowed to percolate through the sorbent using gravity without allowing the sorbent to dry out between steps.

Following the passage of the methanol, DI water, and 0.1 M phosphate buffer (pH 6) through the SPE columns, each diluted sample (that is, calibrator, control, and case item) was loaded on to an individually marked SPE tube and allowed to pass through the sorbent using gravitational flow. The columns were then washed with sequential additions of 3 mL of DI water, 1 mL of 0.1 M acetic acid, and 3 mL of methanol. The SPE columns were then dried by applying a vacuum to the SPE manifold at 15 in. of mercury pressure using an electric vacuum pump.

The analytes were eluted from the SPE columns by the addition of 3 mL of a 78:20:2 methylene chloride–isopropanol–ammonium hydroxide solution. This solution was prepared daily by adding 2 mL of concentrated ammonium hydroxide solution to 20 mL of isopropanol and mixing well. Then, 78 mL of methylene chloride was added to this solution and

Endurance.

It's at the core of every great achievement – and it's the most critical requirement in liquid chromatography pump sealing. We knew this when we set out to design the new, ultra-reliable Enduris™ UHPLC Seal.

Tested to over 2 million cycles at a barrier-breaking 22,000 psi, Enduris combines proven Canted Coil Spring technology, advanced polymer formulation and a precision-engineered backup to deliver pump sealing performance that lasts. And lasts.



ENDURIS[™]
UHPLC SEAL



Marathon testing capabilities, faster throughput, increased sample resolution and fewer service calls. Experience the Enduris performance advantage at www.balseal.com/enduris, or call 800.366.1006




Figure 3: Chromatogram of hair analysis containing amphetamines and synthetic cathinones at the LOQ (0.1 ng/mg). See Table I for compound identification.




Figure 4: Chromatogram of genuine hair demonstrating the presence of amphetamine only (upper trace). Lower trace indicates internal standards. See Table I for compound identification.

the resultant solution was transferred to a clean screw-top bottle for use. A screw-top bottle ensures that the basicity of the solution remains high by eliminating any loss of ammonia from the bottle. The elution solvent was allowed to flow through the SPE sorbent with the aid of gravity and collected in separate glass tubes (75 mm × 12 mm). Glass tubes were chosen because they are standard laboratory materials within this toxicology laboratory.

Following elution, 200 μ L of a solution containing 95% of mobile-phase A and 5% of mobile-phase B was added to

each sample tube. The tubes were vortex mixed for approximately 1 min before the eluates were evaporated to the mobile phase using a gentle stream of nitrogen at 35 $^{\circ}$ C. After further vortex mixing, the samples were transferred to an autosampler vial (2 mL) containing a low-volume insert (250 μ L) and the vial was capped for analysis.

Matrix Effects

Studies of the matrix effects were performed according to a previously published procedure (23). In this part of

the study, aliquots of amphetamine and synthetic cathinones (covering the linear range) were introduced into 200 μ L of a solution containing 95% of mobile-phase A and 5% mobile-phase B. Each elution solution was evaporated to remove the organic component until only the mobile phase remained and then they were analyzed by LC–MS–MS (analysis A). Concurrently, a set of hair samples were subjected to the SPE process; after elution of the analytes from the SPE columns, the elution solvent was spiked with amphetamine, and 200 μ L of mobile phase (95:5 A–B) was added, and then the solution was evaporated to the mobile phase (analysis B). A second set of hair samples was spiked with amphetamine and synthetic cathinones and processed via the SPE method. After elution, 200 μ L of mobile phase was added and after vortex mixing, the solution was evaporated to the mobile phase (analysis C). The data (peak areas) for analyses A, B, and C were collected by the data analysis software. By comparing the peak areas of analysis B with those of A, an assessment of matrix effects was made. The comparison of peak areas for C with B provided data for the recoveries.

Amphetamine and synthetic cathinone solution (concentration: 50 ng/mg) was infused into the tandem mass spectrometer using the on-board syringe pump (controlled by the data analysis software) via a Hamilton syringe (model 1001TLL, 1 mL volume, supplied by Fisher Scientific) at a flow rate of 5 μ L/min. At the same time as the amphetamine solution was flowing into the mass spectrometer, a 10- μ L aliquot of the SPE-extracted hair matrix (samples of hair confirmed to contain no drug material) was injected using the autosampler syringe on the Shimadzu liquid chromatograph. The liquid chromatograph and mass spectrometer were arranged so that samples from the liquid chromatograph were mixed into the flow of amphetamine and synthetic cathinones via a three-port tee section before the total flow entered the mass spectrometer. Any suppression effects on the amphetamine could be monitored at the MRM step for the noted drugs.

Selectivity

When analyzing samples of hair extracts via SPE and LC–MS–MS it is essential to ensure that the interfering effects of

other drug compounds can be eliminated. In this procedure, samples of hair extracts were spiked with 49 drugs at a concentration equivalent to 100 ng/mg of hair sample: (bupropion, lidocaine, methadone, amitriptyline, nortriptyline, thioridazine, trazodone, mesoridazine, meperidine, diphenhydramine, phenyltoloxamine, imipramine, desipramine, benztropine, trimethoprim, diltiazem, haloperidol, strychnine, morphine, codeine, 6-acetylmorphine, oxycodone, oxymorphone, hydrocodone, noroxycodone, hydromorphone, diazepam, nordiazepam, oxazepam, temazepam, alprazolam, α -hydroxyalprazolam, lorazepam, triazolam, α -hydroxytriazolam, flunitrazepam, 7-amino-flunitrazepam, chlordiazepoxide, midazolam, α -hydroxymidazolam, flurazepam, desalkyl-flurazepam, cocaine, ecgonine methyl ester, ecgonine ethyl ester, benzoylecgonine, cocaethylene, clonazepam, and 7-amino-clonazepam) and extracted according to the SPE method. It was observed that the interfering effect of these compounds was not found to be significant.

Results and Discussion

Recovery

The mean recovery of amphetamine and synthetic cathinones from drug-free hair samples was determined to be 95% ($\pm 2\%$). This is an excellent indicator for the efficiency of the extraction procedure of amphetamines and synthetic cathinones using hair as a matrix. This procedure was performed twice daily over a period of five days.

Imprecision of Analysis

The results of the analysis of the spiked control samples of hair (2 ng/mg and 8 ng/mg, respectively) are shown in Table II. Analysis of the control samples was performed at the same time as the calibration curves were constructed — that is, over a period of five days. Control samples were prepared by adding the amphetamine and synthetic cathinone solution to the hair sample (10 mg) in the digestion mixture and treating as per the test samples.

Intraday and interday variation for the analysis of amphetamines and synthetic cathinones was found to be less than 7%

and less than 10%, respectively. Ion suppression studies revealed that suppression of the monitored ions was less than 5%. This method was found to be linear ($r^2 > 0.995$) over the dynamic range 0.1–10 ng/mg.

Limit of Detection and Limit of Quantification

The *limit of detection* (LOD) of a particular method can be defined as the level at which the signal-to-noise ratio for the particular analyte is greater than or equal to 3:1. The *limit of quantification* (LOQ) for the method is the level at which the signal-to-noise ratio for a particular analyte is greater than or equal to 10:1. In this study, LOD values were determined empirically by analyzing extracted samples of drug-free hair fortified with amphetamines and synthetic cathinones by LC-MS-MS according to the SPE method. This was performed until the lowest level at which each of the respective analytes just failed the signal-to-noise ratio of 3:1. This was observed as 0.05 ng/mg. In terms of LOQ, samples of drug-free hair samples were spiked with amphetamines and synthetic cathinones at



SUPELCO
Solutions within.™

BIOshell™ Fused-Core® Columns for Peptides and Proteins

Analyze proteins and peptides by RPLC with high speed and efficiency

BIOshell
U/HPLC Columns

Experience optimum protein recovery and sub-2 μm porous particle efficiency with half of the backpressure.

To learn more about achieving enhanced performance in large molecule applications, visit sigma-aldrich.com/BIOshell

Explore the solutions within



©2014 Sigma-Aldrich Co. LLC. All rights reserved. SAFIC, SIGMA-ALDRICH and SUPELCO are trademarks of Sigma-Aldrich Co. LLC, registered in the US and other countries. Solutions within and BIOshell are trademarks of Sigma-Aldrich Co. LLC. Fused-Core is a registered trademark of Advanced Materials Technology, Inc.

12070

SIGMA-ALDRICH®

concentrations below 10 ng/mg and extracted according to the SPE procedure until the analytes just failed a signal-to-noise ratio of 10:1; this value was found to be 0.1 ng/mg. This is better than the recommendations by the Society of Hair Testing (SoHT) of 0.2 ng/mg (24). Representative chromatograms at LOQ and genuine hair sample are shown in Figures 3 and 4.

Solid-Phase Extraction

Because hair is a solid matrix, the sample requires digestion to produce a liquid that is able to flow through the SPE column. As the digest is basic, buffering to pH 6 permits both efficient flow and optimal sorbing of the drugs onto the SPE sorbent. In using a mixed-mode (C8 and strong cation-exchange chemistries), the sample can be cleaned up via aqueous acid and methanol washes, which leaves the drugs in a much cleaner state than when they were originally applied to the SPE column. This effect is noted in the low matrix effects and ion suppression values.

In this new methodology, the SPE eluates are not evaporated to dryness as in

typical SPE procedures, but are evaporated to a solution of mobile phase. Without the addition of methanolic hydrochloric acid to reduce their volatility, these drugs are known to be lost during this evaporation step. Using the mobile phase as a keeper solvent permits the volatile amphetamines and synthetic cathinones to be kept in solution and results in high recovery values by eliminating the loss during the evaporation.

Tandem MS

In this methodology, LC-MS-MS has been successfully applied to the extraction and analysis of amphetamine and synthetic cathinones rather than GC-MS in which a multistep derivatization procedure (that is, reaction with a fluoroacyl reagent such as heptafluorobutyric anhydride), evaporation, and reconstitution in a volatile solvent is required not only to quantify, but also to confirm the identity of the amphetamine or synthetic cathinone. By using LC-MS-MS with specific MRM values, amphetamine and synthetic cathinones can be targeted, confirmed, and quantified in hair samples without the use

of derivatization. Coupling this procedure with a quick LC method offers analysts the ability to determine concentrations of the drug within a short turnaround time.

Concentrations of Amphetamines and Synthetic Cathinones in Genuine Hair Sample

It was observed that the donor hair sample contained only amphetamine (confirmed by a prescription of Adderall). The concentration was found to be 1.2 ng/mg. The synthetic cathinones (butylone, ethylone, flephedrone, mephedrone, methylone, methedrone, methcathinone [4-MEC], methylenedioxypropylvalerone [MDPV], and pyravalerone) were not observed in the sample and neither were methamphetamine, MDA, or MDMA. Although studies were not performed on different colored hair samples, the recoveries of the drugs on donated hair samples (gray) were notably high.

Conclusion

The use of hair samples is recognized as an alternative specimen of choice for the forensic and clinical analysis of drugs. Analysts

Editors' Series:

Prospects and Pitfalls in Metabolomics: From Study Designs to Interpreting Data in Biomedical Applications

LIVE WEBCAST: Wednesday, March 19, 2014 8:00 am PDT/ 10:00 am CDT/ 11:00 am EDT

Register Free at www.chromatographyonline.com/prospects

EVENT OVERVIEW:

In this web seminar, Professor Oliver Fiehn of UC Davis will explain the importance of different steps in the metabolomic workflow and how to avoid the most dramatic pitfalls along the way. Key topics covered will include defining a suitable hypothesis, designing a study, preparing biological specimen for comprehensive chemical analysis, choosing the most suitable analytical platforms, analyzing raw data including identification of metabolites, using quality controls and validations, and interpreting the data on the background of biomedical questions. Examples will be given from published work using blood plasma in human cohort studies and animal models.

Key Learning Objectives:

- How to choose define a suitable hypothesis for a metabolomics study and design a study
- How to prepare biological specimens for comprehensive chemical analysis
- What analytical platforms to use
- How to analyze raw data
- How to use quality controls and validations



Presenter

Prof. Oliver Fiehn, PhD
Director, West Coast
Metabolomics Center
UC Davis Genome Center



Moderator

Laura Bush
Editorial Director,
LCGC

Who Should Attend:

- Researchers working in metabolomics

Sponsored by



Presented by



For questions, contact Kristen Moore at kmoore@advanstar.com

around the world are going to be asked to test it on a routine basis. With this in mind, this SPE and LC-MS-MS procedure will offer these facilities the ability to perform a quick and efficient analysis of amphetamine and synthetic cathinones in small samples of hair. The novel use of the mobile phase as the keeper solvent serves to improve efficiency and maintain high recovery values.

References

- (1) M. Moawad, M. Khoo, S. Chiang, S. Lee, and J.R. Nennell, *J. AOAC Int.* **93**, 116–122 (2010).
- (2) R.C. Baselt, *Disposition of Toxic Drugs and Chemicals in Man, 7th Edition* (Biomedical Publications, Foster City, California), pp. 66–68 (2004).
- (3) J.F. Wyman, E.S. Lavins, D. Englehart, E.J. Armstrong, K.D. Snell, P.D. Boggs, S.M. Taylor, R.N. Norris, and F.P. Miller, *J. Anal. Toxicol.* **37**, 182–185 (2013).
- (4) B.M. Cawrse, B. Levine, R.A. Jufer, D.A. Fowler, S.P. Vorce, A.J. Dickson, and J.M. Holler, *J. Anal. Toxicol.* **36**, 434–439 (2012).
- (5) K.J. Lusthof, R. Oosting, A. Maes, M. Verschraagen, A. Dijkhuizen, and A.G.A. Sprong, *Forensic Sci. Int.* **206**, e93–95 (2011).
- (6) B.L. Murray, C.M. Murphy, and M.C. Beuhler, *J. Med. Toxicol.* **8**, 69–75 (2012).
- (7) W.A. Baumgartner, V.A. Hill, and W.H. Bland, *J. Forensic Sci.* **34**, 1433–1453 (1989).
- (8) A.H. Beckett and M. Rowland, *J. Pharm. Pharmacol.* **17**, 628–639 (1965).
- (9) P.S. Sever, J. Caldwell, L.G. Dring, and R.T. Williams, *Eur. J. Clin. Pharm.* **6**, 177–80 (1973).
- (10) C.L. Winek, W.W. Wahba, C.L. Winek, Jr., and T.W. Balzer, *Forensic Sci. Int.* **122**, 107–123 (2001).
- (11) Y.H. Caplan and B.A. Goldberger, *J. Anal. Toxicol.* **25**, 396–399 (2001).
- (12) V. Cirimele, P. Kintz, and P. Mangin, *J. Leg. Med.* **108**, 265–267 (1996).
- (13) J.L. Villamour, A.M. Bermejo, P. Fernandez, and M.L. Tabernerero, *J. Anal. Toxicol.* **29**, 135–139 (2005).
- (14) M. Uhl, *Forensic Sci. Int.* **84**, 281–294 (1997).
- (15) C. Jurado, M.P. Gimenez, T. Soriano, M. Menendez, and M. Repetto, *J. Anal. Toxicol.* **24**, 11–16 (2000).
- (16) D. Hensley and J.T. Cody, *J. Anal. Toxicol.* **23**, 518–523 (1999).
- (17) E. Lendoiro, O. Quintela, A. de Castro, A. Cruz, M. Lopez-Rivadulla, and M. Concheiro, *Forensic Sci. Int.* **217**, 207–215 (2012).
- (18) A. de Castro, E. Lendoiro, O. Quintela, M. Concheiro, M. Lopez-Rivadulla, and A. Cruz, *Forensic Toxicol.* **30**, 193–198 (2012).

- (19) *Solid Phase Applications Manual* (UCT Inc., Bristol, Pennsylvania, 2010), pp. 13–15.
- (20) C. Girod and C. Staub, *Forensic Sci. Int.* **107**, 261–271 (2000).
- (21) Y. Gaillard and G. Pepin, *Forensic Sci. Int.* **762**, 251–267 (1997).
- (22) K.A. Hadidi, T. Al-Nsour, and S. Abu-Ragheib, *Forensic Sci. Int.* **135**, 129–136 (2003).
- (23) B.K. Matuszewski, M. Constanzer, and C.M. Chavez-Eng, *Anal. Chem.* **75**, 3019–3030 (2003).
- (24) G.A.A. Cooper, R. Kronstrand, and P. Kintz, *Forensic Sci. Int.* **218**, 20–24 (2012).

Albert A. Elian and Kerrie T. Donovan are with the Massachusetts State Police Crime Laboratory in Sudbury, Massachusetts. **Jeffery Hackett** is with UCT Inc., in Bristol, Pennsylvania. Direct correspondence to: jhackett@unitedchem.com ■

For more information on this topic, please visit www.chromatographyonline.com

EXPERIENCE DIFFERENT

“A Distribution Partner Focused on Your Success”

- ◆ A PARTNER IN SCIENCE
- ◆ CERTIFICATIONS: THAT ASSIST IN DIVERSITY SPEND REQUIREMENTS
- ◆ CAN CO-EXIST WITH PRIMES IN KEY ACCOUNTS
- ◆ 40+ YEARS OF PRACTICAL LAB EXPERIENCE
- ◆ HIGH PROFICIENCY IN SUPPLY CHAIN
- ◆ EFFECTIVE CHANNEL MANAGEMENT SOLUTIONS
- ◆ EXTENSIVE AND GROWING PRODUCT PORTFOLIO
- ◆ FLEXIBLE B2B E-COMMERCE SOLUTIONS
- ◆ PERSONALIZED/CUSTOM PROGRAMS



Bench Top Instruments, Chromatography Supplies, Chemicals, Labware & Life Science Consumables

SBA, SDB, NMSDC, WBE & MBE



PITTCON
CONFERENCE & EXPO 2014

McCormick Place, Chicago, IL
March 3-6, 2014 • Come Visit us at Booth #3916

Contact Us at:

800-343-6015

www.NetaScientific.com

sales@netascientific.com



Diverse Company for Scientific Supplies

Application of Pyrolysis–Gas Chromatography–Mass Spectrometry for the Identification of Polymeric Materials



The analytical pyrolysis technique hyphenated to gas chromatography–mass spectrometry (GC–MS) has extended the range of possible tools for the characterization of synthetic polymers and copolymers. Pyrolysis involves thermal fragmentation of the analytical sample at temperatures of 500–1400 °C. In the presence of an inert gas, reproducible decomposition products characteristic for the original polymer or copolymer sample are formed. The pyrolysis products are chromatographically separated using a fused-silica capillary column and are subsequently identified by interpretation of the obtained mass spectra or by using mass spectra libraries. The analytical technique eliminates the need for pretreatment by performing analyses directly on the solid or liquid polymer sample. In this article, application examples of analytical pyrolysis hyphenated to GC–MS for the identification of different polymeric materials in the plastic and automotive industry, dentistry, and occupational safety are demonstrated. For the first time, results of identification of commercial light-curing dental filling material and a car wrapping foil by pyrolysis–GC–MS are presented.

Peter Kusch*, Gerd Knupp*, Wolfgang Fink*, Dorothee Schroeder-Obst*, Volker Obst†, and Johannes Steinhaus*

*Hochschule Bonn-Rhein-Sieg, University of Applied Sciences, Department of Applied Natural Sciences, Rheinbach, Germany;

†Volker Obst, Dr. Obst Technische Werkstoffe GmbH, Rheinbach, Germany. Direct correspondence to: Peter.Kusch@h-brs.de

Structural analysis and the study of degradation properties are important to understand and improve performance characteristics of synthetic polymers and copolymers in many industrial applications. Traditional analytical techniques used for characterization of polymers and copolymers such as thermal analysis and Fourier transform infrared (FT–IR) spectroscopy have limitations or are not sufficiently sensitive (1). Pyrolysis techniques hyphenated to gas chromatography–mass spectrometry (GC–MS) have extended the range of possible tools for the characterization of synthetic polymers and copolymers. Under controlled conditions, at elevated temperatures (500–1400 °C) in the presence of an inert gas, reproducible decomposition


products characteristic for the original polymer or copolymer sample are formed. The pyrolysis products are chromatographically separated using a fused-silica capillary column and subsequently identified by interpretation of the obtained mass spectra or by using mass spectra libraries (such as the National Institute of Standards and Technology [NIST] or Wiley). Pyrolysis methods eliminate the need for pretreatment by performing analyses directly on the solid polymer or copolymer sample (1). (Please note that this article was presented at the XVII European Conference on Analytical Chemistry, which was held in Warsaw, Poland, on August 25–29, 2013).

Most of the thermal degradation results from free radical reactions



No pre- or post-labeling required

LC-MS analysis for 55 amino acids in 10 minutes






Figure 1: Schematic view of the furnace pyrolyzer used in this study.

initiated by bond breaking and depends on the relative strengths of the bonds that hold the molecules together. A large molecule will break apart and rearrange in a characteristic way (2–4). If the energy transfer to the sample is controlled by temperature, heating rate, and time, the fragmentation pattern is reproducible and characteristic for the original polymer or copolymer. Another sample of the same composition, heated at the same rate to the same temperature for the same period of time, will produce the same decomposition products. Therefore, the essential requirements of the apparatus in analytical pyrolysis are reproducibility of the final pyrolysis temperature, rapid temperature rise, and accurate temperature control. Depending on the heating mechanism, pyrolysis systems have been classified into two groups: the continuous-mode pyrolyzer (furnace pyrolyzer) and pulse-mode pyrolyzer (flash pyrolyzer, such as the heated filament, Curie-point, and laser pyrolyzer). The pyrolysis unit is directly connected to the injector port of a gas chromatograph. A flow of an inert carrier gas, such as helium, flushes the pyrolyzates into the fused-silica




Figure 2: Pyrolysis–GC–MS chromatogram of plastic particles from industrial filter fins at 700 °C obtained with apparatus 1. Fused-silica GC capillary column: 60 m × 0.25 mm, 0.25- μm d_f Elite-5ms. GC conditions: programmed column temperature: 60 °C for 1 min, then 60–100 °C at 2.5 °C/min and then 100–280 °C at 10 °C/min (20 min hold at 280 °C); split-splitless injector temperature: 250 °C; split flow: 50 cm^3/min ; helium programmed pressure: 70 kPa for 1 min, then 70–110 kPa at 1 kPa/min (hold at 110 kPa to the end of analysis). For peak identification, see Table I.

Table I: Pyrolysis products and identified materials in plastic particles from industrial filter fins

Retention Time t_R (min)	Pyrolysis Product at 700 °C	Matching Factor	Identified Material
5.44	Propylene	820	Polypropylene glycol
5.58	1-Butene/1,3-butadiene	840	Styrene–butadiene rubber (SBR)
5.86	Acetone	850	Polypropylene glycol
6.13	Pentadiene	885	SBR
7.43	Benzene	954	SBR
9.72	Toluene	863	SBR
10.51	Cyclopentanone	933	Poly(hexamethylene adipamide) (nylon 6-6)
11.81	2-Cyclopenten-1-one	906	Poly(hexamethylene adipamide) (nylon 6-6)
14.39	Styrene	851	SBR
28.15	4-Isopropylphenol	944	Polycarbonate or bisphenol A epoxy resin
40.40	N-Phenyl-1-naphthalen-amine	948	Antioxidant

capillary column. Figure 1 shows the schematic view of the furnace pyrolyzer used in our investigation. The detection technique of the separated compounds is typically MS, but other GC detectors have also been used depending on the intentions of the analysis (1,4).

The applications of analytical pyrolysis–GC–MS range from research and development of new materials, quality control, characterization and competi-

tor product evaluation, medicine, biology and biotechnology, geology, air-space, and environmental analysis to forensic purposes or conservation and restoration of cultural heritage. These applications cover analysis and identification of polymers, copolymers, and additives in components of automobiles, tires, packaging materials, textile fibers, coatings, half-finished products for electronics, paints or varnishes,




Figure 3: Pyrolysis–GC–MS chromatogram of a car wrapping material at 600 °C obtained with apparatus 1. Fused-silica GC capillary column 60 m × 0.25 mm, 0.25- μ m d_f , Elite-5ms. GC conditions: programmed column temperature: 60 °C for 1 min, then 60–100 °C at 2.5 °C/min and then 100–280 °C at 10 °C/min (20 min hold at 280 °C); split–splitless injector temperature: 250 °C; split flow: 50 cm³/min; helium programmed pressure: 70 kPa for 1 min, then 70–110 kPa at 1 kPa/min (hold at 110 kPa to the end of analysis). For peak identification, see Table II.




Figure 4: Pyrolysis–GC–MS chromatogram of polyaramid fibers at 900 °C obtained with apparatus 2. Fused-silica GC capillary column: 59 m × 0.25 mm, 0.25- μ m d_f , DB-5ms. GC conditions: programmed column temperature: 75 °C for 1 min, then 75–280 °C at 7 °C/min (hold to the end of analysis); programmed pressure of helium carrier gas: 122.2 kPa for 1 min, then 122.2–212.9 kPa at 7 kPa/min (hold at 212.9 kPa to the end of analysis). For peak identification, see Table III.

lacquers, leather, paper or wood products, food, pharmaceuticals, surfactants, and fragrances.

Our earlier publications (1,5–12) presented the analysis and identification of degradation products of commercially available synthetic polymers and copolymers by using analytical pyrolysis hyphenated to gas chromatography with flame ionization detection (GC–FID) and GC–MS. In this work, new examples of applications of this analytical technique for the identification of different polymeric materials are demonstrated.

Experimental Samples

Plastic particles from industrial filter fins, a car wrapping foil, unknown fibers, and commercial light-curing dental filling material were used in the investigation.

Instrumentation and Analytical Conditions

Approximately 100–200 μ g of solid sample was cut out with a scalpel and inserted without any further preparation into the bore of the pyrolysis solids-injector and then placed with

the plunger on the quartz wool of the quartz tube of the furnace pyrolyzer Pyrojector II (SGE Analytical Science). Three spots on each sample were analyzed in duplicate. The pyrolyzer was operated at a constant temperature of 550, 600, 700, or 900 °C. The pressure of helium carrier gas at the inlet to the furnace was 95 kPa.

Pyrolysis–GC–MS measurements were made using two apparatus. In the first apparatus (1), the pyrolyzer was connected to a Trace 2000 gas chromatograph (ThermoQuest, CE Instruments) with a quadrupole mass spectrometer Voyager (ThermoQuest, Finnigan, MassLab Group) operated in electron ionization (EI) mode. A 60 m × 0.25 mm, 0.25- μ m Elite-5ms fused-silica GC capillary column (PerkinElmer Instruments) was used. The GC conditions were as follows: programmed column temperature: 60 °C for 1 min, then 60–100 °C at 2.5 °C/min, 100–280 °C at 10 °C/min (20-min hold at 280 °C). The temperature of the split–splitless injector was 250 °C and the split flow was 50 cm³/min. Helium, grade 5.0 (Westfalen AG), was used as a carrier gas. The helium programmed pressure was 70 kPa for 1 min, then 70–110 kPa at 1 kPa/min (hold at 110 kPa to the end of analysis) was used. The transfer line temperature was 280 °C. The MS EI ion source temperature was kept at 250 °C. The ionization occurred with a kinetic energy of the impacting electrons of 70 eV. The current emission of the rhenium filament was 150 μ A. The MS detector voltage was 350 V. Mass spectra and reconstructed chromatograms (total ion current [TIC]) were obtained by automatic scanning in the mass range m/z 35–450 u. Pyrolysis–GC–MS data were processed with the Xcalibur software (ThermoQuest) and the NIST 05 mass spectra library.

In the second apparatus (2), the pyrolyzer was connected to a 7890A gas chromatograph with a series 5975C quadrupole mass spectrometer (Agilent Technologies Inc.) operated in EI mode. A 59 m × 0.25 mm, 0.25- μ m d_f DB-5ms fused-silica GC capillary column (J&W Scientific) was used. Helium, grade 5.0 (Westfalen AG), was used as a carrier gas. The GC conditions were as follows: programmed

Table II: Pyrolysis products and identified materials in car wrapping foil			
Retention time t_R (min)	Pyrolysis product at 700 °C	Matching factor	Identified material
5.49	Hydrogen chloride	945	Poly(vinyl chloride) (PVC)
5.58	Methyl chloride	800	PVC
5.64	1-Butene	938	PVC
6.01	1,3-Pentadiene	921	PVC
6.99	Tetrahydrofuran	769	Solvent
7.29	1,4-Cyclohexadiene	923	PVC
7.48	Benzene	945	PVC
8.20	Methyl methacrylate	750	Poly(methyl methacrylate) (PMMA)
9.78	Toluene	904	PVC
10.29	2-Ethyl-1-hexene	890	Bis(2-ethylhexyl) phthalate (plasticizer)
10.46	Cyclopentanone	912	Poly(hexamethylene adipamide) (nylon 6-6)
10.53	1-Octene	907	PVC
14.47	Styrene	934	PVC
17.33	2-Ethylhexanal	856	Bis(2-ethylhexyl) phthalate (plasticizer)
20.57	2-Ethyl-1-hexanol	916	Bis(2-ethylhexyl) phthalate (plasticizer)
21.07	<i>o</i> -Methylstyrene	888	PVC
21.43	Indene	870	PVC
22.99	<i>p</i> - <i>tert</i> -Butyltoluene	856	2,6-Bis-(1,1-dimethylethyl)-4-methylphenol (BHT) (antioxidant) (?)
25.93	Naphthalene	920	PVC
28.46	2-Methylnaphthalene	862	PVC
28.63	Phthalic anhydride	906	Bis(2-ethylhexyl) phthalate (plasticizer)
28.79	1-Methylnaphthalene	875	PVC
36.69	3,3-Diphenylacrylonitrile	937	Adhesive layer

column temperature: 60 °C for 1 min, then 60–280 °C at 7 °C/min (hold at 280 °C to the end of analysis); programmed helium pressure: 122.2 kPa for 1 min, then 122.2–212.9 kPa at 7 kPa/min (hold at 212.9 kPa to the end of analysis). Second set of GC conditions: programmed column temperature: 75 °C for 1 min, then 75–280 °C at 7 °C/min (hold at 280 °C to the end of analysis); programmed helium pressure: 122.2 kPa for 1 min, then 122.2–212.9 kPa at 7 kPa/min (hold at 212.9 kPa to the end of analysis).

The temperature of the split-splitless injector was 250 °C and the split ratio was 50:1. The transfer line temperature was 280 °C. The MS EI ion source temperature was kept at 230 °C. The ionization occurred with a kinetic energy of the impacting electrons of 70 eV. The quadrupole temperature was 150 °C. Mass spectra and reconstructed chromatograms (total ion current) were obtained by automatic scanning in the mass range m/z 35–750 u. Pyrolysis–GC–MS data were processed with the ChemStation software (Agilent Technologies) and the NIST 05 mass spectra library.

Results and Discussion

Pyrolysis–GC–MS of Plastic Particles from Industrial Filter Fins

A sample of plastic particles from industrial filter fins was pyrolyzed at 700 °C to identify its composition. Figure 2 shows the obtained pyrolysis–GC–MS chromatogram of the sample. Based on the decomposition products summarized in Table I, the plastic particles were identified as a mixture of poly(hexamethylene adipamide) (nylon 6-6) and polypropylene glycol with a small amount of styrene–butadiene–rubber (SBR). The peaks of propylene and acetone indicate the presence of polypropylene glycol. The main decomposition product of nylon 6-6 is cyclopentanone (retention time $[t_R] = 10.51$ min). Other peaks in Figure 2, like butene/1,3-butadiene ($t_R = 5.58$ min), benzene ($t_R = 7.43$ min), toluene ($t_R = 9.72$ min), and styrene ($t_R = 14.39$ min) are typical pyrolysis products of SBR (1,2,5,6,18). The small peak of 4-isopropylphenol ($t_R = 28.15$ min) may be a clue to the presence of polycarbonate or bisphenol A epoxy resin (5,6). All of the pyrolysis products and the materials




Figure 5: Chemical structure of poly(*p*-phenylene terephthalamide) (polyaramid).

identified from pyrolysis products in filter fins are summarized in Table I.

Pyrolysis–GC–MS of a Car Wrapping Foil

The next object of identification was a car wrapping foil pyrolyzed at 600 °C. Figure 3 shows the obtained pyrolysis–GC–MS chromatogram of the car wrapping foil. Based on the decomposition products summarized in Table II, the plastic material was identified as a mixture of flexible poly(vinyl chloride) (PVC) with bis(2-ethylhexyl) phthalate (BEHP) plasticizer and poly(hexamethylene adipamide) (nylon 6-6). The chromatogram in Figure 3 shows the typical pyrolysis products of PVC, like hydrogen chloride ($t_R = 5.49$ min), benzene ($t_R = 7.48$ min), and naphthalene ($t_R = 25.93$ min). This is the result of the formation of double bonds by the elimination of hydrogen chloride from the poly(vinyl chloride) macromolecules, followed by the breaking of the carbon chain with or without cyclization reaction (2). The detected cyclopentanone ($t_R = 10.46$ min) is generally known as characteristic pyrolysis product of nylon 6-6 (2,3,6). Methyl methacrylate ($t_R = 8.20$ min) identified in pyrolyzate is formed from poly(methyl methacrylate) (6) and most likely comes from the adhesive film. Thus, the identified 3,3-diphenylacrylonitrile ($t_R = 36.69$ min) may be from the adhesive layer of the foil.

The thermal decomposition of the plasticizer bis(2-ethylhexyl) phthalate identified in car wrapping foil leads to the formation at 600 °C of 2-ethyl-1-hexene ($t_R = 10.29$ min), 2-ethylhexanal ($t_R = 17.33$ min), 2-ethyl-1-hexanol ($t_R = 20.57$ min), and phthalic anhydride ($t_R = 28.63$ min) (1,7). In the car wrapping material, the rest of the tetrahydrofuran solvent ($t_R = 6.99$ min) was also detected. Table II shows the identified ingredients of the pyrolyzed car wrapping foil.

Identification of Unknown Plastic Fibers

A sample of unknown plastic fibers was pyrolyzed at 700 °C and 900 °C, respectively, to identify its composition. Figure 4 shows the pyrolysis–GC–MS chromatogram of the sample pyrolyzed at 900 °C. Based on the decomposition products summarized in Table III, the

fibers were identified as polyaramid [poly(*p*-phenylene terephthalamide)] (Figure 5). The main identified degradation products of polyaramid at 900 °C are benzene ($t_R = 7.56$), aniline ($t_R = 11.61$ min), and benzonitrile ($t_R = 11.81$ min). Currently, polyaramid fibers have only been characterized in a few publications using thermal analysis (thermogravimetry, derivative thermogravimetry, and differential thermal analysis), infrared spectroscopy techniques (13–15), and pyrolysis–GC–MS (2,16–19).

Polyaramid fibers are a class of heat-resistant, strong synthetic fibers. They are used in aerospace and military applications for ballistic-rated body armor, fabric, ballistic composites, and fire fighters protective clothing as well as in bicycle tires and as an asbestos substitute.

Identification of Commercial Light-Curing Dental Filling Material

A number of dental filling materials are presently available for tooth restorations.



Automate Your Food/Flavor/Fragrance Analysis

Introducing EST Analytical's FLEX GC Autosampler

Now you're in control.

Sample Preparation and Sample Introduction

- Liquid, Headspace and SPME capability
- Robust and easy-to-use automation
- Customizable
- On-board sample prep routines
- Cost-effective



Discover a new level of laboratory automation designed to improve your everyday analytical life. *That's purposeful innovation... that's*

EST 
analytical

503 Commercial Drive, Fairfield, OH 45014
800.283.3510 | estanalytical.com

Table III: Pyrolysis products of polyaramid fibers

Peak Number	Retention Time t_R (min)	Pyrolysis Product at 900 °C	Matching Factor
1	6.52	Carbon dioxide	957
2	6.70	Acrylonitrile	896
3	7.56	Benzene	968
4	8.50	Toluene	918
5	10.20	Styrene	924
6	11.30	Isocyanatobenzene	947
7	11.61	Aniline	959
8	11.81	Benzonitrile	967
9	16.98	1,2-Benzodinitrile	959
10	17.02	1,4-Benzodiamine	905
11	19.60	Biphenyl	957
12	26.19	Acridine	931
13	26.64	1,1'-Biphenyl-4-amine	960
14	28.31	Carbazole	928
15	30.01	N-Phenylbenzamide	916




Figure 6: Pyrolysis–GC–MS chromatogram of commercial light-curing dental filling material at 550 °C obtained with apparatus 2. Fused-silica GC capillary column: 59 m \times 0.25 mm, 0.25- μ m d_f DB-5ms. GC conditions: programmed column temperature: 60 °C for 1 min, then 60–280 °C at 7 °C/min (hold at 280 °C to the end of analysis); programmed helium pressure: 122.2 kPa for 1 min, then 122.2–212.9 kPa at 7 kPa/min (hold at 212.9 kPa to the end of analysis). For peak identification, see Table IV.

The four main groups of these materials, which dentists have used for about 35 years, are the conventional glass-ionomer cements, resin-based composites, resin-modified glass-ionomer cements, and polyacid-modified resinous composites (20). Light-curing glass-ionomer cements contain polyacrylic acid, chemically or photo-curing monomers (multifunctional methacrylates, like triethylene glycol dimethacrylate or 2-hydroxyethyl methacrylate), an ion-leaching glass, and additives (initiators,

inhibitors, stabilizers, and others) (20). Resin-modified glass-ionomer cements are now widely used in dentistry as direct filling materials, liners, bases, luting cements, and fissure sealants (21). These materials mainly consist of polymer matrix and glass-ionomer parts. The polymer matrix is based on a monomer system and different multifunctional methacrylates with additives (21). Methacrylic monomers, such as bisphenol A glycidyl methacrylate (Bis-GMA), urethane dimethacrylate

(UDMA), triethylene glycol dimethacrylate (TEGDMA), and 2-hydroxyethyl methacrylate (HEMA), are the main components of resin-based dental filling materials. The presence of additives such as initiators, activators, inhibitors, and plasticizers in uncured dental material mixture is necessary (21).

Figure 6 shows the total ion current pyrolysis–GC–MS chromatogram of commercial light-curing dental filling material pyrolyzed at 550 °C. The pyrolysis products identified by using mass spectra library NIST 05 are summarized in Table IV. The carbon dioxide (t_R = 6.85 min) identified in pyrolyzate is formed from polyacrylic acid (2,18). The identified substances 2-hydroxyethyl methacrylate (HEMA) (t_R = 13.65 min), ethylene glycol dimethacrylate (EGDMA) (t_R = 19.48 min), and triethylene glycol dimethacrylate (TEDMA) (t_R = 28.72 min) are known as standard composites of dental filling materials (1). Other compounds in Table IV, such as bisphenol A (t_R = 33.10 min) or bisphenol A diglycidyl ether (t_R = 42.42 min), are probably formed by thermal degradation of bisphenol A diglycidyl mono- or dimethacrylates. The presence of the additives, such as the antioxidant 2,6-bis(1,1-dimethylethyl)-4-methylphenol (BHT) (t_R = 23.17 min) or the UV-absorber drometrizol (t_R = 31.95 min) was also confirmed. The triphenylantimony (t_R = 34.55 min) identified in pyrolyzate is used as catalyst in the UV-induced polymerization (1).

Conclusion

Analytical pyrolysis–GC–MS has been proven as a valuable technique for the analysis and identification of organic polymeric materials in the plastic and rubber industry. For the first time pyrolysis–GC–MS was used for the identification of commercial light-curing dental filling material and for the identification of a car wrapping foil. This technique allows the direct analysis of very small sample amounts (5–200 μ g) without the need for time-consuming sample preparation.

References

- (1) P. Kusch, in *Advanced Gas Chromatography – Progress in Agricultural, Biomedical and*

Table IV: Pyrolysis products of commercial light-curing dental filling material

Peak Number	Retention Time t_R (min)	Pyrolysis Product at 550 °C	Matching Factor
1	6.85	Carbon dioxide	999
2	9.62	Methacrylic acid	936
3	12.96	Phenol	909
4	13.65	2-Hydroxyethyl methacrylate (HEMA)	918
5	19.40	4-Isopropenylphenol	928
6	19.48	Ethylene glycol dimethacrylate (EGDMA)	915
7	23.00	Not identified	
8	23.17	2,6-Bis-(1,1-dimethylethyl)-4-methylphenol (BHT)	923
9	23.65	Not identified	
10	23.89	Not identified	
11	28.72	Triethylene glycol dimethacrylate (TEDMA)	958
12	31.95	Drometrizol (Tinuvin-P)	938
13	33.10	4,4'-Dihydroxy-2,2-diphenylpropane (bisphenol A)	924
14	34.55	Triphenylantimony	911
15	35.25	Tetraethylene glycol dimethacrylate	791
16	36.98	Not identified	
17	42.42	Bisphenol A diglycidyl ether (BADGE)	839

Industrial Applications, M.A. Mohd, Ed. (InTech, Rijeka, Croatia, 2012), pp. 343–362.

- S.C. Moldoveanu, *Analytical Pyrolysis of Synthetic Organic Polymers* (Elsevier, Amsterdam, The Netherlands, 2005).
- T.P. Wampler, *Applied Pyrolysis Handbook, 2nd Ed.* (CRC Press, Boca Raton, Florida, 2007).
- K.L. Sobehi, M. Baron, and J. Gonzales-Rodrigues, *J. Chromatogr. A* **1186**, 51–66 (2008).
- P. Kusch, *Chem. Anal.* (Warsaw) **41**, 241–252 (1996).
- P. Kusch, G. Knupp, and A. Morrisson, in *Horizons in Polymer Research*, R.K. Bregg, Ed. (Nova Science Publishers, New York, New York, 2005), pp. 141–191.
- P. Kusch, *LCGC North Am.* **31**(3), 248–254 (2013).
- P. Kusch, V. Obst, D. Schroeder-Obst, W. Fink, G. Knupp, and J. Steinhaus, *Eng. Fail. Anal.* **35**, 114–124 (2013).
- P. Kusch and G. Knupp, *Nachr. Chem.* **57**, 682–685 (2009).
- P. Kusch, V. Obst, D. Schroeder-Obst, G. Knupp, and W. Fink, *LCGC AdS*, July/August, 5–11 (2008).
- P. Kusch and G. Knupp, *LCGC AdS*, June, 28–34 (2007).
- P. Kusch, W. Fink, D. Schroeder-Obst, and V. Obst, *Aluminium* (Isernhagen, Germany) **84**(4), 76–79 (2008).
- S. Villar-Rodil, A. Martínez-Alonso, and J.M.D. Tascón, *J. Anal. Appl. Pyrol.* **58–59**, 105–115 (2001).
- F. Suárez-García, A. Martínez-Alonso, and J.M.D. Tascón, *Carbon* **42**, 1419–1426 (2004).

- S. Villar-Rodil, A. Martínez-Alonso, and J.M.D. Tascón, *J. Therm. Anal. Calorim.* **79**, 529–532 (2005).
- J.R. Brown and A.J. Power, *Polym. Degrad. Stab.* **4**(5), 379–392 (1982).
- H.-R. Schulten, B. Plage, H. Ohtani, and S. Tsuge, *Angew. Makromol. Chem.* **155**, 1–20 (1987).
- S. Tsuge, H. Ohtani, and C. Watanabe, *Pyrolysis-GC/MS Data Book of Synthetic Polymers* (Elsevier, Amsterdam, The Netherlands, 2011).
- Zs. Czégény and M. Blazsó, *J. Anal. Appl. Pyrol.* **58–59**, 95–104 (2001).
- R. Rogalewicz, A. Voelkel, and I. Kownacki, *J. Environ. Monit.* **8**, 377–383 (2006).
- R. Rogalewicz, K. Batko, and A. Voelkel, *J. Environ. Monit.* **8**, 750–758 (2006).

Peter Kusch, Gerd Knupp, Wolfgang Fink, Dorothee Schroeder-Obst, and Johannes Steinhaus are with Hochschule Bonn-Rhein-Sieg, University of Applied Sciences in the Department of Applied Natural Sciences, in Rheinbach, Germany. **Volker Obst** is with Dr. Obst Technische Werkstoffe GmbH, in Rheinbach, Germany.

Direct correspondence to: Peter.Kusch@h-brs.de ■

For more information on this topic, please visit www.chromatographyonline.com

SAVE THE DATE

35th Annual

Minnesota Chromatography Forum Spring Symposium

May 7th & 8th, 2014
Earle Brown Heritage Center
Minneapolis, MN

Advanced Registration Ends April 30, 2014.


For more information, visit: www.minnchrom.com

- Short Courses (May 6 & 7)
- Vendor Seminars (May 7)
- Exhibition (May 7 & 8)
- Technical Presentations (May 8)
- Keynote Speaker (May 8)

Professor Jonathan Sweedler

Dept of Chemistry, University of Illinois.

Topic: "MS-Based Characterization Of The Brain: From Individual Cells To Cellular Release"



The Minnesota Chromatography Forum is a scientific society that was founded in 1978 to maintain and promote education, discussion and exchange of information with respect to all fields of chromatography.

PRODUCTS & RESOURCES

MS system

The ionKey/MS system from Waters is designed for applications such as bioanalyses, pharmacokinetics studies, food safety analyses, and environmental analyses. According to the company, the system is intended for use with the company's Acquity UPLC M-Class system and Xevo TQS mass spectrometer.

Waters Corporation,
Milford, MA.

www.waters.com/ionkey



SEC system

Agilent's 1260 Infinity Multi-Detector Bio SEC Solution system is designed for the analysis of proteins and biotherapeutics. According to the company, the instrument includes a photodiode-array detector and a light-scattering detector, bioinert flowpaths, and intuitive software to determine accurate molecular weight, size, and aggregate presence with high reproducibility.

Agilent Technologies,
Santa Clara, CA.

www.agilent.com/chem/infinity-bio-sec




HPLC columns

Diamond-based Flare HPLC columns from Diamond Analytics are designed to expand the existing range of analytical capabilities in separation. According to the company, the columns are stable at pH 11.3 and 120 °C and can be regenerated for repeat use.

Diamond Analytics,
Orem, UT.

www.diamond-analytics.com



HPLC and UHPLC columns product bulletin

A product bulletin from Advanced Materials Technology contains information about the company's BioClass reversed-phase HPLC and UHPLC columns. According to the company, the eight-page bulletin includes nine figures that illustrate the benefits of the columns. The bulletin also includes product specifications and part numbers.

Advanced Materials Technology, Inc.,

Wilmington, DE.

www.advanced-materials-tech.com



GC autosampler

The Flex Series GC autosampler from EST Analytical is designed for both end users and OEM applications. According to the company, the autosampler has liquid injection capability with an upgrade path to headspace or SPME analysis.

EST Analytical,
Fairfield, OH.

www.estanalytical.com



GC system

The TRACE 1300 gas chromatography system from Thermo Scientific is designed to conserve helium, permit use of multiple detectors simultaneously, automate sampling gas workflows, and perform flame photometric detection.

According to the company, new options include a helium-saver auxiliary oven, a dedicated flame photometric detector, and a gas sampling valve module. **Thermo Fisher Scientific,** Santa Clara, CA.

www.thermoscientific.com/trace1300



Thermal desorption autosampler

The 7500 TD autosampler from CDS Analytical is designed as an extension of the company's 7400 purge-and-trap autosampler. According to the company, the autosampler has the ability to change between TD tubes for air and VOA vials for soil and water analysis.

CDS Analytical,
Oxford, PA.

www.cdsanalytical.com



GC valves

Analytical Flow Products GC valves from Norgren are designed to prevent cross-contamination and excessive wear in challenging applications. According to the company, a purge groove design improves data quality, extends instrument life, increases environmental safety, and decreases downtime.

Norgren,

Littleton, CO.

www.norgren.com



LC | GC's CHROMacademy

powered by crawfordscientific

Infrared Spectroscopy Training

2 New eLearning Courses

Available Now at CHROMacademy

www.chromacademy.com/infrared



Gas generator

The Corona Nitrogen 1010 gas generator from Peak Scientific is designed to supply pure nitrogen to Thermo Fisher Scientific's Corona Veo charged aerosol detector. According to the company, the generator has been fully tested and approved.

Peak Scientific,

Billerica, MA.

www.peakscientific.com




LC columns

Restek's Raptor ARC-18 superficially porous particle LC columns are designed for use with LC-MS-MS systems. According to the company, the columns can be used in low pH (≤ 2.0) conditions and provide consistent retention, peak shape, and response for charged bases, neutral acids, and small polar compounds.

Restek Corporation,

Bellefonte, PA.

www.restek.com



HPLC columns

Selectra HPLC columns from UCT include a pentafluorophenylpropyl phase, a polyaromatic phase, and a C18 phase. According to the company, the columns are produced using high purity, type B spherical silica, the material is fully endcapped, and columns are available with either 3- μm or 5- μm particles.

UCT, LLC,

Bristol, PA.

www.unitedchem.com



GPC instrument

Tosoh Bioscience's Eco-SEC high-temperature GPC system is designed to be used at temperatures ranging from 40 °C to 220 °C. According to the company, the dual-pump system includes a dual-flow refractive index detector, a column oven, system temperature control, an autosampler, a sample prep station, and control software.

Tosoh Bioscience, LLC, King of Prussia, PA.

www.tosohbioscience.com




SEC-MALS detector

The mDAWN multiangle light scattering (MALS) detector from Wyatt Technology is designed to be coupled to any UHPLC system to determine absolute molecular weights and sizes of polymers, peptides, and proteins or other biopolymers directly, without using column calibration or reference standards. According to the company, the instrument's flow cell volume is less than 10 mL.

Wyatt Technology Corp., Santa Barbara, CA.

www.wyatt.com



Photodiode-array detector

Shimadzu's SPD-M30A photodiode-array detector is designed for a variety of HPLC and UHPLC conditions and reportedly can be used for a range of analyses without replacing its capillary cell. According to the company, the detector's capillary cell allows the peak from the principal component and a 0.005% infinitesimal peak to be quantified simultaneously.

Shimadzu Scientific Instruments,

Columbia, MD.

www.ssi.shimadzu.com




LC columns

ACE UltraCore solid-core LC columns from Advanced Chromatography Technologies are designed to provide a low column back pressure and are available in SuperC18 and Super-PhenylHexyl bonding. According to the company, both phases feature proprietary encapsulated bonding technology for peak shape and phase stability across a pH range of 1.5 to 11.0.

Advanced Chromatography Technologies Ltd,


Aberdeen, Scotland. www.ace-hplc.com




Report development and generation software

Bruker Dash Reporting software, from Bruker Chemical and Applied Markets is designed to provide customized reporting that centers on Dash Designer, a purpose-built standalone application that allows users to position and closely format report elements, and preview reports with relevant data. According to the company, individual elements can be sorted, filtered, resized, and formatted with common editing operations and advanced functions.

Bruker Corporation, Fremont, CA. www.bruker.com





Multi-Antibiotic Residue Detection in Food:

An Improved Method for Screening and Confirmation Testing, in Accordance with EU Commission Decision 2002/657/EC

ON-DEMAND WEBCAST

Register Free at: www.chromatographyonline.com/detection

EVENT OVERVIEW:

In the complex field of testing food for veterinary drug residues, innovative known/unknown screening techniques using high-resolution accurate mass (HRAM) LC-MS offer unique opportunities and provide improved cost efficiency. Compared to conventional screening methods such as inhibitory tests or other biotests, HRAM can cover a large set of targeted substances, including transformation or metabolic products, with high selectivity and sensitivity.

Tune into this web seminar to learn about a newly developed HRAM-LC-MS/MS method for screening, quantitation, and confirmation of more than 100 antibiotic, anthelmintic, and antiparasitic residues in a single run. The method has been fully validated in accordance with 2002/657/EC guidelines and can be applied for analysis of veterinary drug residues in samples of animal origin, such as meat or fish.

Key Learning Objectives:

- Discover how high-resolution accurate mass spectrometry helps to overcome traditional challenges of MS/MS based methods
- Learn how to improve your lab throughput and save costs by performing screening, quantitation, and confirmation of residues in a single run
- Learn how to improve confidence in your data by eliminating false positives and negatives

For questions, contact Kristen Moore at kmoore@advanstar.com



PRESENTER

Nelli Jochim, PhD
Mass Spectrometry Specialist
Eurofins WEJ Contaminants GmbH



MODERATOR

Laura Bush
Editorial Director
LC/GC

Who Should Attend:

- Routine laboratory managers and technical specialists working in residue analysis
- Analytical professionals from regulatory and official control authorities
- Analytical chemists interested in new developments in accurate mass spectrometry

Sponsored by

Thermo
SCIENTIFIC

Presented by

LC|GC

Food safety analysis virtual conference

EMD Millipore's Lab Solutions Virtual Conference on Food Safety Analysis reportedly includes three webinars: "ISO Standardization in Food Microbiology and Quality of Media," "The Importance of Water Quality in Food Analyses," and "Monolith Chromolith Columns – an Ideal Tool for the Analysis of Food Samples with Complex 'Dirty' Matrices." Registration for on-demand viewing is available at <https://engage.vevent.com/index.jsp?eid=2484&seid=425>

EMD Millipore, Billerica, MA. www.millipore.com



Custom HPLC columns

Custom HPLC columns from Hamilton are designed to allow any of the company's stationary phases to be packed into almost any hardware format. According to the company, the HPLC stationary phases are available in varying particle sizes and in most column hardware dimensions.

Hamilton, Company, Reno, NV.
www.hamiltoncompany.com



Quick-connect system

The Opti-Lynx II system from Optimize Technologies is designed as a combination of quick-connect holders with a selection of packed-bed cartridges. According to the company, accessing and changing the insert takes a quarter turn and the connection is rated up to 6000 psi.

Optimize Technologies, Oregon City, OR.
www.optimize-tech.com



Sorbents

Supelco's Supel QuE Z-Sep, Z-Sep/C18, and Z-Sep+ sorbents are designed to provide robust LC-MS and GC-MS methods for a variety of analytes in difficult matrices. According to the company, the Z-Sep family of sorbents is available in 2-mL and 12-mL tub formats for QuEChERS.

Supelco/Sigma-Aldrich, Bellefonte, PA.
www.sigma-aldrich.com/zsep



Sample cleanup workstation

The Freestyle workstation from Pickering Laboratories is designed for automated sample cleanup work flow. The instrument is based on a suspended rack design, with an XY robot arm for liquid handling. The workstation reportedly is able to handle multiple flask shapes with volumes ranging from 1 mL to 1 L, and the instrument's software enables users to program multiple sample parameters and to prepare graphical reports and audit logs. SPE, GPC, and evaporation and solvent exchange modules are available.

Pickering Laboratories, Inc., Mountain View, CA. www.pickeringlabs.com



Hydrogen gas generators

High-pressure hydrogen generators from Parker Hannifin are designed to produce up to 1300 cc/min of 99.99999+% pure hydrogen and eliminate cylinders of hydrogen fuel gases and helium carrier gases. According to the company, a single generator can support up to 20 instruments with fuel and carrier gas.

Parker Hannifin, Haverhill, MA.
www.parker.com/fns/balstonlabgas-generators



Capillary tubing

Polymicro's flexible fused-silica capillary tubing is designed with an outer diameter of 1/32 in. According to the company, the tubing mates with existing 1/32-in. fittings and is available in a range of internal diameters from 50 µm to 500 µm.

Polymicro Technologies, Phoenix, AZ.
www.polymicro.com



Peptide mapping columns

AdvanceBio peptide mapping columns from Neta Scientific are designed for resolution and identification of amino acid modifications in primary structure. According to the company, the columns feature a 120-Å pore size with superficially porous 2.7-µm particles.

Neta Scientific, Hainesport, NJ.
www.netascientific.com



The Use of Solid Core Technology in the Pharmaceutical Environment

LIVE WEBCAST: Tuesday, March 25, 8:00 am PDT/ 11:00 am EDT/ 3:00 pm GMT

Register free at www.chromatographyonline.com/core_technology

EVENT OVERVIEW:

Since the introduction of solid core technology, there has been great interest within the separation science field to look at ways in which this innovative technology can be employed in a variety of different industries. One such area is the pharmaceutical arena, where there is a broad range of sample types as well as requirements throughout the process of developing new chemical entities. This presentation will look at how the technology of solid core can be readily adapted to cope with the various challenges associated with the pharmaceutical sector, looking at various sample matrices and various molecular entities, from small molecules to large biomolecules.

The presentation will seek to advise users on appropriate column selection, with an insight into how varying the solid core to porous layer can allow the user to optimize separation performance by reducing extra band broadening. Data will also be presented that demonstrates how this technology is inherently more robust than standard fully porous systems when analyzing biological extracts, routinely used in DMPK departments, resulting in longer column lifetimes.

Presented by

LC|GC

Sponsored by

Thermo
SCIENTIFIC

Key Learning Objectives:

- How to choose an appropriate solid core column dependent on the application
- How solid core morphology can be optimized for small and large molecule analysis
- How solid core technology can extend column lifetimes

Who Should Attend:

- Anybody involved in analyzing pharmaceutical compounds

Presenters



Dafydd Milton
Product Manager,
LC & LC-MS Columns
Chromatography
Consumables
Thermo Fisher Scientific



Tony Edge
R&D Principal
Chromatography
Consumables
Thermo Fisher Scientific



Moderator

Laura Bush
Editorial Director
LCGC

For questions, contact Kristen Moore at kmoore@advanstar.com

AD INDEX

ADVERTISER	PAGE	ADVERTISER	PAGE
Advanced Materials Technology	198	OI Analytical	199
Agilent Technologies	CV2, 224	Optimize Technologies	201
Bal Seal Engineering Company	205	Pace Analytical Services	217
Bruker Daltonics	165	Parker Hannifin Corporation	CV3
CDS Analytical, Inc.	168	Peak Scientific Instruments Ltd.	185
CASSS	192, 196	Pickering Laboratories	168
Chem Service, Inc.	164	Polymicro Technologies LLC	181
Detector Engineering & Technology	164	Regis Technologies, Inc.	183
EMD Millipore	171	Restek Corporation	169
EST Analytical	215	Shimadzu Scientific Instruments	191
Hamilton Company	163	Spark Holland	187
Hichrom	193	Supelco, Inc.	CV TIP, 170, 207
Imtakt USA	211	Thermo Fisher Scientific	189, 195, 221, 223
LECO Corporation	208	Tosoh Bioscience	177
Metrohm USA	CV4	UCT, LLC	179
Nacalai USA, Inc.	197	Waters Corporation	167
Neta Scientific	209	Wyatt Technology Corporation	173
Norgren	175		

Sample Preparation Method Development for **Unconventional Matrices**

LIVE WEBCAST: Wednesday, March 12, 1:00 pm EDT

Register Free at www.chromatographyonline.com/matrices

EVENT OVERVIEW:

Are you faced with developing or improving your sample preparation methods for unusual matrices? If so, this seminar will provide a simple, straightforward approach to modifications that will increase compound recovery from unique matrices focusing on QuEChERS and Solid Phase Extraction methodology.

QuEChERS (Quick, Easy, Cheap, Effective, Rugged and Safe) was developed to be a simple, effective, and affordable way to extract pesticide residues from a wide range of fruit and vegetable matrices. The benefits of QuEChERS are well documented, making this

technique very popular and easily executed in the laboratory with common equipment such as pipettes, vortexer, and centrifuge. Because of the benefits QuEChERS provides as a sample preparation technique, it is rapidly expanding to other matrices and compound classes; it's just not for fruits and vegetables anymore. Therefore, simple modifications might be required to the basic procedure to optimize compound recovery from these unique matrices.

Sometimes a simple technique like QuEChERS does not offer the degree of matrix removal or LLOQ required. In that case, a solid phase extraction (SPE) method is implemented to provide ultimate cleanliness and recovery. Method development for SPE is thought to be an arduous task. However, with a simple step-by-step approach, optimization of a SPE method is easily obtained resulting in the desired results.

- Evaluating the components that make a sample matrix complex
- Optimizing a QuEChERS method for unique matrices and compound classes
- Optimizing SPE parameters in order to achieve highest recovery and cleanliness

Key Learning Objectives:

- To be able to evaluate a unique sample matrix and understand the parameters that will need to be adjusted, evaluated
- How to systematically adjust parameters easily
- Getting the most out of your sample preparation technique: QuEChERS and SPE

Who Should Attend:

- Chemists and lab managers performing GC, LC, or mass spectrometry for food, forensic, environmental, chemical, pharmaceutical, cosmetic, personal care, and consumer goods testing, who need to improve chromatographic results and increase lab productivity by streamlining their sample preparation processes.



Presenter

Joan Stevens, Ph.D.
Senior Applications Scientist,
Sample Preparation Products
Agilent Technologies



Moderator

Laura Bush
Editorial Director,
LCGC

Sponsored by



Agilent Technologies

Presented by

LCGC



Products • Services • Equipment
Training • Career Opportunities

FROM

LCGC
north america




QbD (QUALITY BY DESIGN):
VOLS I AND II COMBO **\$299.90**

QbD (Quality by Design): A systematic approach to product and process design and development, compiles the best content from BioPharm International and Pharmaceutical Technology to provide valuable insight into the topic and assist you in making the business case for QbD based on the criteria decision makers need to evaluate initiatives and related technology.

ELEMENTS OF BIOPHARMACEUTICAL
PRODUCTION SERIES, FOURTH EDITION

**Newly
Updated**



Led by Dr. Anurag Rathore. For anyone involved in or planning to start process development, characterization and/or validation activities.

\$129.95

BioPharm INTERNATIONAL

LCGC north america

Pharmaceutical
Technology

save on these and other educational resources at
www.industrymatter.com | 800.598.6008



industrymatter

comprehensive information for professionals

H P L C INSTRUMENTS & SUPPLIES

SONNTEK RX

Your Prescription for
Quality Research Lamps!

201-236-9300

www.sonntek.com

sonntek@aol.com

Call Tod McCloskey
at 800-225-4569 ext. 2739
to place your ad!
tmccloskey@advanstar.com



THE ESSENTIALS

Excerpts from LCGC's professional development platform, CHROMacademy.com

Mixed-Mode HPLC Separations: What, Why, and How

Although reversed-phase columns (such as C18) are used in a wide range of applications, they often fail to retain highly polar (or charged) analytes, and offer limited selectivity.

Hydrophilic-interaction chromatography (HILIC) is often used to analyze polar analytes, but it encounters challenges such as increased use of organic solvents, poor solubility of analytes in highly organic media, retention affects created by the sample matrix, and a limitation to the extent to which the retention of hydrophobic analytes can be controlled.

Ion-exchange chromatography can be used to retain charged molecules, but it fails to retain neutral analytes and can show poor selectivity for analytes of the same charge.

A current estimate of the usage of these techniques is reversed phase ~60+%; HILIC or normal phase ~18%; ion exchange ~5%.

The use of mixed-mode stationary phases can offer a solution to some of the problems outlined above and their properties and applications can be summarized as follows:

- Mixed-mode chromatography is a method that uses more than one separation mode; mainly reversed-phase combined with ion-exchange interactions, which allows the retention and separation of both polar and nonpolar analytes in a single analysis.
- The biggest benefit of this approach is that selectivity can be optimized by adjusting mobile-phase ionic strength, pH, or organic solvent. As a result, the selectivity can be finely tuned for the separation of compounds with widely different physico-chemical properties. For example, drug molecules and their counterions may be separated in a single analysis.
- Mixed-mode chromatography requires no ion-pairing agents in the mobile phase

for separating highly hydrophilic charged analytes, which simplifies the mobile phase and is compatible with mass spectrometry (MS).

A general schematic for the various approaches to producing mixed-mode stationary phases is shown in Figure 1.

The mixed-beads approach in Figure 1 is a blend of two different stationary phases into a single column: reversed phase (in black) and ion exchange (in red). Issues with this approach include concerns regarding the batch-to-batch reproducibility and homogeneity of the resulting high performance liquid chromatography (HPLC) columns.

The mixed-ligand approach to mixed-mode stationary phase production was a second-generation product, but again, ensuring the absolute reproducibility and homogeneity of the surface is problematic.

The embedded and tipped ligands are a third-generation approach to mixed-mode phases and produce homogeneous and reproducible mixed-mode phases. The position of the ion-exchange group within the bonded phase ligand dictates the relative strength of the hydrophobic and electrostatic interactions and the embedded ligand can be described as an ion-exchange–modified reversed-phase column and the tipped ligand as a reversed-phase–modified ion-exchange column.

The zwitterionic ligand is a fourth-generation approach to mixed-mode separations in which both anionic and cationic functional groups are included in each ligand and the strength and position of the ion-exchange

moieties can be altered to affect the selectivity of each ligand towards various analytes.

In mixed-mode separations, the hydrophobic moieties of analyte molecules are retained by interaction with the hydrophobic parts of the bonded ligand, and the charged groups are either retained or repelled by the ion-exchange groups on the stationary phase surface. There are several ways to control retention and selectivity with mixed-mode phases, as one might imagine.

The amount of organic modifier can be altered to affect changes in retention and selectivity of the separation, primarily by changing the hydrophobicity of the eluent, as would be typical for a reversed-phase separation.

Eluent pH can also be used to alter selectivity and retention, and the pH will affect not only the degree of ionization of the analyte molecule, but also the extent of the surface charge if weak anion- or cation-exchange ligands are used. The ionic strength of the eluent can be used to alter the separation characteristics, and ionic-strength gradients are often used to increase the elution strength of the eluent with respect to the electrostatic interactions between analytes and the stationary phase surface.

The nature of mixed-mode stationary phases offer several possibilities to “tune” the selectivity of a separation and as such offer an important tool in the analytical armory for anyone dealing with highly polar analytes or analytes whose physico-chemical properties differ widely.




Figure 1: Schematic representation of the various approaches to producing mixed-mode stationary phases.

More Online:

Get the full tutorial at

www.CHROMacademy.com/Essentials
(free until April 20).



Together, we can conserve nonrenewable helium by converting GCs to hydrogen.

Helium is a non-renewable resource. There is a global shortage of helium that has caused prices to increase significantly. Many GC users are considering switching to hydrogen as a carrier gas.

Hydrogen provides many advantages over helium including, higher resolution, shorter run times, longer column life, and cost savings. Many instrument manufacturers now provide resources to make switching carrier gas easier.

A hydrogen generator is a safe, convenient, and inexpensive solution to storing hydrogen cylinders. Learn more at solutions.parker.com/conservhelium.

aerospace
climate control
electromechanical
filtration
fluid & gas handling
hydraulics
pneumatics
process control
sealing & shielding

Find your IC Greatness

Because when switching to Metrohm you'll get a great system, backed by even greater support.

Guaranteed!



Your methods will run – right out of the box!



You're trained to be the IC pro in your lab



You'll enjoy expert care & 72 hr on site service!



Get a free T-shirt (no purchase necessary)!

Complete details at:
www.ic-changeisgood.com/guaranteed



866-METROHM (638-7646)

